

spirit of excellence

SOLF



Real-time fluorescence imaging with the ENDOCAM[®] Logic 4K camera platform

System green

Real-time fluorescence imaging with the ENDOCAM Logic 4K camera platform

System green from Richard Wolf is a visualization solution for real-time ICG/NIR fluorescence imaging which meets the highest standards of endoscopic imaging.

ICG/NIR fluorescence imaging has become one of the most valuable technologies in the field of minimally invasive surgery for applications including intraoperative assessment of visceral perfusion (e.g., anastomosis), visualization of biliary structures in fluorescence cholangiography, and lymph node mapping (e.g., visualization of sentinel lymph nodes during lymph node dissection).

High-resolution image quality meets ground-breaking visualization

System green, an extension of the ENDOCAM Logic 4K camera platform, enables brilliant fluorescence imaging in real time.

- Simultaneous ICG/NIR and white light endoscopy with brilliant image quality
- Clear visualization of visceral perfusion
- Superior exposure of biliary anatomy in fluorescent cholangiography
- Precise visualization of lymphatic system





System green can also be used for conventional, high-resolution "white light" endoscopy. With only a touch of a button, users can flip between white light and fluorescence imaging.




ENDOLIGHT LEDgreen

Powerful illumination for clear visualization

ENDOLIGHT LEDgreen is an innovative solution for use in both standard white light endoscopy and ICG/NIR fluorescence imaging.

The use of advanced LED technology eliminates the need to use a laser diode to visualize the ICG/NIR information. The LED illumination produces a homogeneous, consistently bright endoscopic image allowing the user to cover two applications with one system.

- Easy operation—select white light or ICG/NIR illumination with a simple press of a button
- dialog Technology—automatic adjustment of light intensity in conjunction with ENDOCAM Logic HD & 4K camera controllers
- Reduced service and operating costs as a result of extended lifespan of latest LED technology (up to 15,000 hours)



The absorption spectrum of ICG (indocyanine green) is located in the near-infrared range (NIR). ICG absorbs the light in the range of 600 nm to 900 nm and emits the fluorescence between 800 nm up to 850 nm. Thus, **System green** allows visualization and detection of the fluorescence intraoperatively. core



ß SOLF N Profile greenICG -LAP/THORAX WOLF greenICG greenICG pure • **ENDOCAM** Logic 4K ß 1 V • ENDOLIGHT .ED green



In greenICG mode, the conventional white light endoscopic image is simultaneously combined with the ICG/NIR image information in real time, at the touch of a button. This allows the tissue to be differentiated with a high level of detail and contrast.



System green

Real-time fluorescence imaging with the ENDOCAM Logic 4K camera platform



ENDOCAM Logic 4K Camera Controller

The ENDOCAM Logic 4K camera controller enables endoscopic imaging in native 4K UHD resolution with ICG/NIR visualization via two unique Special Imaging Modes.

- Special Imaging Modes:
 - greenICG (white light image with ICG overlay)
 - greenICG pure (grayscale image with ICG overlay)
- High-resolution 4K UHD endoscopic image
- Existing camera controllers can easily be upgraded to the latest software enabling ICG/NIR capability

System green is designed to meet the needs of the busy Operating Room. Each component guarantees efficient, flexible, and future-proof use across multidisciplinary laparoscopic applications.

ENDOCAM Logic Camera Head green

- For use with ENDOCAM Logic HD and 4K camera controllers
- Also suitable for use in white light endoscopy
- Digital image processing and signal transmission with dedicated powerful NIR sensor for high fluorescence sensitivity
- Lightweight and ergonomic

rodic areau

 Integrated lens (f=24 mm) optimized for use in laparoscopy



spirit of excellence

PANOVIEW ULTRA Telescopes

With improved fiber optics and a superior lens system, the PANOVIEW ULTRA Telescopes deliver a stunning level of brightness and clarity.

- Optimized for 4K UHD and ICG/NIR applications
- Increased Number of Light Fibers deliver more light to the surgical site than traditional fiber optic technologies
- Low-Dispersion Glass Rod Lenses for improved peripheral sharpness
- Superior fiber optics for natural color visualization in both white light and System green mode—greenICG


Fusion_Fiber Light Cables

Engineered with fused-glass technology, which eliminates space between individual light fibers, fusion light cables transmit approximately 20% more light than standard light cables.

- Higher temperature resistance
- Anti-kink protection
- Autoclavable for cost-effective sterilization
- Increased light transmission compared to standard light cables



The greenICG pure mode enables specific, high-contrast visualization as well as differentiation of tissue. The grayscale visualization enables detailed examination of key structures (e.g., lymphatic vessels and lymph nodes).



System green



Camera Controller

ENDOCAM Logic 4K Camera Controller Bundle 552530115 Includes:

ENDOCAM Logic 4K camera controller (5525301), remote control (5525401), power cable (N710171), 8 GB USB stick (56540028), HDMI/ HDMI cable, 3.0 m (103847)

Camera Head

ENDOCAM Logic Camera Head green with integrated lens (f=24mm) 85525812

Telescopes

PANOVIEW ULTRA Teles	cope,
10 mm, 305 mm WL	
0°	
30°	8934462
50°	. 8934463
PANOVIEW ULTRA Teles	cope,

10 mm, 440 mm WL	it made
0°	89344416
30°	89344426
50°	89344436

(N710171) Foot control rockerswitch for ENDOLIGHT LEDgreen.....2030.105

ENDOLIGHT LEDgreen Bundle

ENDOLIGHT LEDgreen light source

(5165002), patch cable, RJ45 SFTP,

0.5 m (72325378), power cable

Light Source

Includes:

Fusion_Fiber Light Cable

adapter projector-side (8095.07), adapter endoscope-side (809509) and fusion fiber light cable (80655030), 3.0 m

4K Medical Monitors

27" Sony (Richard Wolf-Branded) Monitor.....LMD-X2705N 27" Sony (Richard Wolf-Branded) Monitor Bundle

.....LMD-X2705N_BNDL Includes:

27" Sony (Richard Wolf-Branded) Monitor (LMD-X2705N), monitor cover (N800406), thumb screws (N900692) and power supply (AC-300MD)

32" Sony (Richard Wolf-Branded) Monitor.....LMD-X3200N 32" Sony (Richard Wolf-Branded) Monitor Bundle

.....LMD-X3200N_BNDL Includes:

32" Sony (Richard Wolf-Branded) Monitor (LMD-X3200N), monitor cover (53700031), thumb screws (N900692) and power supply (AC-300MD)

31" Foreseeson (FSN) Monitor FS-P3102D 55" Sony (Richard Wolf-Branded) Monitor. LMD-X550NB

Monitor Accessories

Power Supply.	AC-300MD*
Thumb Screws (Pack	(of 4)
	N900692
27" Richard Wolf-Bra	inded
Monitor Cover	N800406
32" Richard Wolf-Bra	Inded
Monitor Cover	53700031
55" Richard Wolf-Bra	inded
Monitor Cover	53700055

Richard Wolf Medical Instruments 353 Corporate Woods Parkway Vernon Hills, IL 60061-3110

Sony Medical Monitors

*For use only with 27", 32" and 55"

(800) 323-WOLF (9653) info@richardwolfusa.com

www.richard-wolf.com/en-us/



spirit of excellence