



DR. WEIGERT

Systematic Hygiene



spirit of excellence

Richard Wolf
Medical Instruments Corp.
353 Corporate Woods Parkway
Vernon Hills, IL 60061 - U.S.A.
Phone: 847-913-1113

Medical



neodisher®
MultiZym

With triple enzyme power!

Multi-enzymatic detergent for reprocessing thermostable and thermolabile instruments in immersion and ultrasonic baths.

neodisher® MultiZym – Active Enzymes for Manual Instrument Reprocessing.

neodisher® MultiZym

Multi-enzymatic detergent with protease, lipase and amylase

- Intended for use for the manual cleaning of thermo-stable and thermolabile instruments, incl. MIS and micro instruments, dental instruments, flexible endoscopes, anaesthetic equipment, containers and other medical utensils in immersion baths and ultrasonic baths
- **neodisher® MultiZym** is not a hazardous product according to the Appendix D, OSHA Hazard Communication Standard 29 CFR 1910:1200



1.32 liquid gallon
(5 l) canister
DI-442735




12 x 33.8 fluid ounce
(1 l) bottle
DI-442744

Our formula for success:

- The formula of **neodisher® MultiZym** combines among others three classes of enzyme, thus guaranteeing maximum cleaning power with a synergistic effect for all types of residue: protease breaks down proteins from protein residues, lipase removes fat-containing contaminants and amylase removes polysaccharides, for example from bronchial mucus.

Application recommendation:

Manual cleaning with normal levels of contamination	0.125 oz/gal - 0.375 oz/gal, 59 - 122 °F, 2 - 10 min 1 - 2.5 ml/l, 15 - 50 °C, 2 - 10 min
Manual cleaning with high levels of contamination	0.375 oz/gal - 3.75 oz/gal, 59 - 122 °F, 2 - 10 min 2.5 - 30 ml/l, 15 - 50 °C, 2 - 10 min
Manual cleaning in ultrasonic baths	0.125 oz/gal - 3.75 oz/gal, 59 - 122 °F, 2 - 10 min 1 - 30 ml/l, 15 - 50 °C, 2 - 10 min



Enzymes work in accordance with the 'lock and key' principle: each enzyme cleaves one specific substrate. Like a key, the enzyme's structure 'fits' one particular type of substrate.

Our Competence – Your Benefit:

- first-class, residue-free cleaning performance
- multi-enzymatic formula with protease, lipase and amylase
- free of boric acid and borates, perfume-free with no colourings; instruments are clearly visible in the solution
- practical one-litre bottle with integrated dosing aid
- independently assures excellent decontamination with regard to organic materials, prevents redeposition of protein residues
- effectiveness of active ingredients, high efficiency and exceptional cleaning results even at low concentrations
- very user-friendly with exceptional material compatibility: suitable for surgical stainless steel, optics, standard synthetic materials, anodised aluminium, silicone and materials used in anaesthetic equipment
- easy to use and simple dosing without additional accessories



Manufactured by:
Chemische Fabrik Dr. Weigert GmbH & Co. KG • Mühlenhagen 85 • Germany-20539 Hamburg • Tel.: +49-40-7 89 60-0 • Fax: +49-40-7 89 60-120 • info@drweigert.de • www.drweigert.com
Distributed by:
Richard Wolf Medical Instruments Corp. • 353 Corporate Woods Parkway • Vernon Hills, IL 60061 - U.S.A. • Phone: 847-913-1113

TO ORDER, PLEASE CONTACT:
Richard Wolf Medical Instruments Corp. • 847-913-1113 • richardwolfusa.com
For professional use only.