

Where ergonomic and intuitive  
design come together



**ERAGON** *axial*

An axial handle for a modular system

# ERAGON<sub>axial</sub>

Where ergonomic and intuitive design come together

**ERAGON** axial handles are the most recent addition to the **ERAGON** modular forceps line. Four new axial handle options offer a fatigue-free operational platform and precision with a new level of handling.

**ERAGON** axial combines the benefits of ergonomics, balance, maneuverability and secure locking functionality. Feel for yourself.



## ○ Ergonomic design gains the upper hand

Designed for user-friendly operation, **ERAGON** axial handles incorporate angled contours that permit a relaxed working position in the hands of surgeons. Wide contact surfaces around the eyelets offer competent fingers a fatigue-free operational platform, alleviating pressure points and maximizing comfort.



○ **Turning with ease**

An optional silicone attachment extends the rotation knob toward the user to ensure the right fit for all hand sizes. The attachment allows for easier access and operation and can easily be removed for reprocessing.

**Safety under control**

The intuitive ratchet lock provides sensitive handling and a secure grip.

**Ultimate compatibility**

- **ERAGON axial** is fully compatible with all components in the ERAGON modular series. All sheaths and jaw inserts can be combined as required.
- The **"Click-it"** system provides the same intuitive assembly and disassembly features of other **ERAGON** modular handles.

**Make the most of the multi-faceted and diverse ERAGON modular series.**

# ERAGON<sub>axial</sub>

Where ergonomic and intuitive design come together



ERAGON axial handle  
with ratchet lock  
**83930082**



ERAGON axial handle  
without lock,  
with HF connector  
**83930083**



ERAGON axial handle  
with ratchet lock,  
with HF connector  
**83930084**






ERAGON axial handle  
with ratchet lock,  
not rotatable  
**83930085**



Option:  
Rotation Extension  
ERAGON axial,  
autoclavable  
**8988**

## ERAGON<sub>modular</sub>

All ERAGON modular components  
are compatible with ERAGON axial.

Catalog Number	Description	WL
8394933	Bariatric 	450 mm
8393933	Standard 	330 mm
8392933	Small anatomy 	260 mm

