

**Simple, low-cost alternative
to active smoke evacuation.**



Smoke Evacuation Filter

Safely captures and filters surgical smoke
and potentially hazardous by-products.

Smoke Evacuation Filter

The simple, low-cost alternative to active smoke evacuation.



spirit of excellence

The Richard Wolf Smoke Evacuation Filter offers a safe, efficient, and sterile option for evacuating smoke and odor generated when using HF current during laparoscopic interventions. Surgical smoke may contain harmful particulates and other by-products, presenting a potential health hazard to OR staff.

- Combined active charcoal and ULPA filter captures and filters 99.999% of particles, smoke odors, and other potentially hazardous by-products.*
- Sterile and individually packed.
- Easy connection to trocars with a 3-mm luer lock.
- Significantly reduces odor.
- Helps to improve the visual field in the abdomen.
- Filters particulates and other potentially hazardous by-products as small as 0.027 microns that may be present in the smoke.



Smoke Evacuation Filter

Pack of 10, disposable product, individually packed, sterile

4170504

Specifications subject to change without notice.

1912-07-01-0619USA

*Data on file; available on request.

To arrange a trial, contact your local Richard Wolf representative or call Customer Service at (800) 323-9653.

Richard Wolf Medical Instruments
353 Corporate Woods Parkway
Vernon Hills, IL 60061-3110

(800) 323-WOLF (9653)
info@richardwolfusa.com

www.richardwolfusa.com