


The laparoscopic toolbox.




ERAGON modular

Laparoscopic instrument system
for multidisciplinary use.

ERAGON^{modular}

The laparoscopic toolbox—a system for laparoscopic specialties.

A tray of ERAGONmodular instruments holds versatility and quality. These modular instruments aid a wide range of surgical applications, pairing a surgeon's expertise with precision, reliability, and intuitive handling.



3.5 mm, 5 mm, and 10 mm

The ERAGONmodular line is available in a wide variety of jaw patterns with varying diameters and working lengths for laparoscopic applications. Ergonomic handles complement a surgeon's preferred sheath and jaw pattern configuration.

Ergonomic Handles

ERAGONmodular handles feature wide contact surfaces at sensitive points, optimum operation of the rotation dial, precise grip and manipulation, and a refined ratchet mechanism.

3.5 mm,



Safety First

The ERAGONmodular sheaths are insulated with Halar®, a high-quality layer of insulation. The Halar® insulation is less likely to sustain damage as compared to shrink tubing insulation. The Halar® coating provides a reliable bond with the metal sheath, as well as enhanced protections when using monopolar HF current. The interior positioning of the locking mechanism for the sheath and handle also safeguard against inadvertent mechanical injury to surrounding tissue.


Ease of Use

The "Click-it" lock technology of the ERAGONmodular system enables easy assembly and disassembly by hospital staff.

Overview of ERAGONmodular

- Reusable instrumentation that is cost-effective
- Intuitive assembly
- Ergonomic design
- Halar® insulation
- Surgeon's preference
 - Wide range of jaw patterns with different diameters and working lengths
 - Large selection of handles
- For minimally invasive applications

Validated for sterilization while assembled or disassembled.



5 mm, and 10 mm

ERAGON_{modular}

The way to an individualized laparoscopic tray.

The Toolbox Principle

Customizable sets for minimally invasive procedures. ERAGON_{modular} offers surgeons a wide range of options for creating the ideal laparoscopic instrument set. The toolbox principle streamlines the tray to the surgeon's preference, reducing costs to the hospital and minimizing the amount of instrumentation needed.

Select Handle


Select Working Length

Select Sheath Tube


Select Jaw Insert




∅
3.5 mm



∅
5 mm




∅
10 mm



















Complete Instrument Set and Tray











ERAGON modular mini

ø 3.5 mm, Insulated

		Axial Handle			
					
All jaw inserts and sheaths can also be combined with ERAGONaxial.		With Lock	Without Lock, with HF Connector	With Lock, with HF Connector	With Lock, Not Rotatable
Extension Attachment	Product No.	83930082	83930083	83930084	83930085
	8988	•	•	•	





Jaw Insert	Sheath Tube, ø 3.5 mm, Insulated	Handle		Complete Instrument • Jaw Insert • Sheath Tube • Handle		
						
	WL 240 mm	WL 330 mm	Without Lock, with HF, Mono	With Lock, with HF, Mono		
	Product No.	8390933	8391933	83930073	83930074	Bundle No.
Grasping and dissecting forceps	8390207	•			•	83902077
 Curved left, double-action	8391207		•		•	83912077
Grasping and dissecting forceps	8390221	•			•	83902217
 Long, fine tip, double-action	8391221		•		•	83912217
Grasping and dissecting forceps	8390222	•			•	83902227
 Short, fine tip, double-action	8391222		•		•	83912227
Grasping and dissecting forceps	8390223	•			•	83902237
 Curved, double-action	8391223		•		•	83912237








∅ 3.5 mm, Insulated

Jaw Insert		Sheath Tube, ∅ 3.5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
		WL 240 mm	WL 330 mm	Without Lock, with HF, Mono	With Lock, with HF, Mono	Without Lock, with HF Connector	
							
Product No.		8390933	8391933	83930073	83930074	83930083	Bundle No.
	8390215	•		•			83902157
	8391215		•	•			83912157
Double-action							
	8390224	•		•			83902247
	8391224		•	•			83912247
Fine serrated, curved, double-action							
	8390227	•		•			83902277
	8391227		•	•			83912277
Single-action							
	8390216	•			•		83902167
	8391216		•		•		83912167
	8390217	•			•		83902177
	8391217		•		•		83912177








Complete instruments in the above table are **recommendations only**. They are comprised of the components: jaw insert, sheath, and compatible handle. However, **any combination is possible within the same diameter**. The ERAGONaxial handle can also be selected as an addition or an alternative.

∅ 3.5 mm, Non-Insulated

		Axial Handle	
 All jaw inserts and sheaths can also be combined with ERAGONaxial.			
		With Lock, without HF Connector	With Lock, without HF Connector, Not Rotatable
Extension Attachment	Product No.	83930082	83930085
	8988	•	







Jaw Insert		Sheath Tube, ∅ 3.5 mm, Non-Insulated		Handle	Complete Instrument • Jaw Insert • Sheath Tube • Handle
					
		WL 240 mm	WL 330 mm	With Lock, without HF	
	Product No.	8390934	8391934	83930072	Bundle No.
Grasping and dissecting forceps	8390219	•		•	83902197
 Double-action	8391219		•	•	83912197
Grasping and dissecting forceps	8390220	•		•	83902207
 2 x 4 teeth, double-action	8391220		•	•	83912207
Atraumatic grasping forceps	8390208	•		•	83902087
 Fenestrated, double-action	8391208		•	•	83912087
Atraumatic grasping forceps	8390218	•		•	83902187
 Fenestrated, single-action	8391218		•	•	83912187






∅ 3.5 mm, Non-Insulated

Jaw Insert		Sheath Tube, ∅ 3.5 mm, Non-Insulated		Handle	Complete Instrument • Jaw Insert • Sheath Tube • Handle
					
		WL 240 mm	WL 330 mm	With Lock, without HF	
	Product No.	8390934	8391934	83930072	Bundle No.
Atraumatic grasping forceps 	8390225	•		•	83902257
Double fenestrated, curved, single-action	8391225		•	•	83912257
Atraumatic grasping forceps 	8390226	•		•	83902267
Fenestrated, double-action	8391226		•	•	83912267
Universal grasping forceps 	8390209	•		•	83902097
Double-action	8391209		•	•	83912097
Grasping forceps "Babcock" 	8390210	•		•	83902107
Double-action	8391210		•	•	83912107












Complete instruments in the above table are **recommendations only**. They are comprised of the components: jaw insert, sheath, and compatible handle. However, **any combination is possible within the same diameter**. The ERAGONaxial handle can also be selected as an addition or an alternative.

∅ 5 mm, Insulated

		Axial Handle			
					
All jaw inserts and sheaths can also be combined with ERAGONaxial.		With Lock	Without Lock, with HF Connector	With Lock, with HF Connector	With Lock, Not Rotatable
Extension Attachment	Product No.	83930082	83930083	83930084	83930085
	8988	•	•	•	












Jaw Insert	Sheath Tube, ∅ 5 mm, Insulated			Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle	
	WL 260 mm	WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono		
	Product No.	8392933	8393933	8394933	83930072	83930073	83930074	Bundle No.
Atraumatic grasping forceps	8392.181	•					•	83921817
	8393.181		•				•	83931817
Fenestrated, with cross-cut toothing, double-action	8394.181			•			•	83941817
Plate-shaped grasping forceps								
	8393.182		•				•	83931827
Fenestrated with fine cross-cut toothing, double-action								
Grasping and dissecting forceps								
	8393.184		•				•	83931847
Fenestrated, with fine cross-cut toothing, 26 mm long, single-action	8394.184			•			•	83941847
Grasping and dissecting forceps								
	8393.185		•				•	83931857
Fenestrated, with fine cross-cut toothing, 21 mm long, single-action								

∅ 5 mm, Insulated












Jaw Insert		Sheath Tube, ∅ 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
							
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
	Product No.	8393933	8394933	83930072	83930073	83930074	Bundle No.
Grasping forceps "Babcock"  Grasping surface with fine cross-cut toothing, double-action	8393.186	•		•			83931867
Grasping forceps "Babcock" 	8393.187	•		•			83931877
Large distal grip jaws with fine cross-cut toothing, double-action	8394.187		•	•			83941877
Atraumatic grasping forceps 	8393.195	•			•		83931957
Smooth, double-action	8394.195		•		•		83941957
Grasping and dissecting forceps  Atraumatic, fine cross toothing, double-action	8393210	•				•	83932107
Grasping and dissecting forceps 	8393220	•				•	83932207
Atraumatic, double-action	8394220		•			•	83942207
Grasping forceps 	8393223	•		•			83932237
Atraumatic, fine cross toothing, double-action	8394223		•	•			83942237

Complete instruments in the above table are **recommendations only**. They are comprised of the components: jaw insert, sheath, and compatible handle. However, **any combination is possible within the same diameter**. The ERAGONaxial handle can also be selected as an addition or an alternative.







ø 5 mm, Insulated

Jaw Insert		Sheath Tube, ø 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
							
	Product No.	WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	Bundle No.
Grasping forceps "Babcock" 	8393242	•		•			83932427
Fenestrated, single-action	8394242		•	•			83942427
Grasping forceps 	8393250	•		•			83932507
Atraumatic, fenestrated, fine cross toothing, double-action	8394250		•	•			83942507
Grasping and coagulating forceps 	8393.287	•			•		83932877
Distal cross-cut toothing with cupped jaw, double-action							
Atraumatic grasping forceps 	8393.292	•		•			83932925
Cupped jaw with wavy tooth edge, 26 mm, single-action	8394.292		•	•			839429255
Atraumatic grasping forceps 	8393.293	•		•			83932935
Cupped jaw with wavy tooth edge, 26 mm, double-action	8394.293		•	•			839429355
Universal grasping forceps 	8393.481	•				•	83934817
Pyramid-shaped toothing and cross-cut toothing, 18 mm, double-action							

∅ 5 mm, Insulated












Jaw Insert		Sheath Tube, ∅ 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
							
Product No.		8393933	8394933	83930072	83930073	83930074	Bundle No.
Universal holding forceps  Pyramid-shaped toothing and cross-cut toothing, 26 mm, double-action		8393.482					83934827
Grasping and extraction forceps "Pike-Mouth"  Double-action		8393.483					83934837
Grasping forceps  Double-action, with row of teeth		8393.485					83934857
Grasping forceps  One tooth, fine cross-cut toothing, slimline, double-action		8393.486					83934867
Grasping forceps  Two teeth, fine cross-cut toothing, single-action		8393.487					83934877
Grasping and extraction forceps  Myoma grasping forceps, 2 x 3-claw, 29 mm, single-action		8393.488					83934887

∅ 5 mm, Insulated














Jaw Insert		Sheath Tube, ∅ 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
							
	Product No.	8393933	8394933	83930072	83930073	83930074	Bundle No.
Grasping forceps							
	8393.490	•		•			83934907
Two teeth, wavy grasping surface, double-action	8394.490		•	•			83944907
Claw grasping forceps							
	8393.494	•		•			83934945
Double-action							
Grasping and extraction forceps							
	8393.188	•		•			83931887
Triangular, fenestrated, with fine cross-cut toothing, double-action, with spring jaws							
Grasping and dissecting forceps							
	8393.189	•		•			83931897
Fine cross-cut toothing, fenestrated, double-action, with spring jaws	8394.189		•	•			83941897
Grasping and dissecting forceps "Debaquey"							
	8393.190	•		•			83931907
Atraumatic clip, double-action, jaw section length 39 mm							

∅ 5 mm, Insulated

spirit of excellence











Jaw Insert		Sheath Tube, ∅ 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
							
Product No.		8393933	8394933	83930072	83930073	83930074	Bundle No.
Grasping and dissecting forceps "Debakey"  Atraumatic clip, double-action, jaw section length 48 mm		8393191					83931917
Atraumatic grasping forceps  Double window, curved, single-action		8393229					83932297
Atraumatic intestinal grasping forceps  Single fenestration, double-action		8393249					83932497
Grasping forceps "Debakey"  Curved, jaw length 35 mm, double-action		8393256					83932567
Grasping forceps "Debakey"  Double-action, jaw length 47 mm		8393258					83932587
Grasping and dissecting forceps "Debakey"  Atraumatic clip longitudinal grooves with fine cross-cut toothings, double-action, with spring jaws		8393.291					83932917

ø 5 mm, Insulated

Jaw Insert		Sheath Tube, ø 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
	Product No.	8393933	8394933	83930072	83930073	83930074	Bundle No.
Grasping and dissecting forceps  Fine cross-cut toothing, double-action	8393.183	•			•		83931837
Grasping and dissecting forceps  Right-angled, double-action	8393192	•				•	83931927
Grasping and dissecting forceps  Jaw length 16 mm, double-action	8393193	•		•			83931937
Grasping and dissecting forceps  Jaw length 16 mm, double-action	8394193		•	•			83941937
Grasping and dissecting forceps  Long, jaw length 24 mm, double-action	8393194	•				•	83931947
Grasping and dissecting forceps  Long, jaw length 24 mm, double-action	8394194		•			•	83941947
Grasping and dissecting forceps  Long, jaw length 30 mm, double-action	8393197	•				•	83931977
Grasping and dissecting forceps  Long, jaw length 30 mm, double-action	8394197		•			•	83941977
Grasping forceps  Curved, fine tip, double-action	8393198	•				•	*
Grasping forceps  Curved, fine tip, double-action	8394198		•			•	83941987
Grasping and dissecting forceps  Long, fine tip, double-action	8393199	•				•	83931997
Grasping and dissecting forceps  Long, fine tip, double-action	8394199		•			•	83941997












*No bundle number available. Components need to be quoted separately (jaw insert, sheath, handle).

∅ 5 mm, Insulated











Jaw Insert		Sheath Tube, ∅ 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
							
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
	Product No.	8393933	8394933	83930072	83930073	83930074	Bundle No.
Grasping and dissecting forceps							
	8393200	•				•	83932007
Short, fine tip, double-action	8394200		•			•	83942007
Grasping and dissecting forceps							
	8393230	•		•			83932307
Short rounded jaws, double-action	8394230		•	•			83942307
Grasping and dissecting forceps							
	8393240	•				•	83932407
Curved, double-action	8394240		•			•	83942407
Grasping and dissecting forceps "Maryland Dissector"							
	8393.281	•				•	83932817
Curved, with fine cross-cut toothing, double-action	8394.281		•			•	83942817
Dissecting forceps "Petelin"							
	8393.282	•				•	83932827
30° left angled, with cross-cut toothing, single-action							

Complete instruments in the above table are **recommendations only**. They are comprised of the components: jaw insert, sheath, and compatible handle. However, **any combination is possible within the same diameter**. The ERAGONaxial handle can also be selected as an addition or an alternative.

ø 5 mm, Insulated











Jaw Insert		Sheath Tube, ø 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
							
Product No.		8393933	8394933	83930072	83930073	83930074	Bundle No.
 Grasping and dissecting forceps "Mixer"	8393.283	•				•	83932837
Angled, fine pyramid-shaped toothing, double-action							
 Grasping and dissecting forceps "Dolphin"	8393.285	•			•		83932857
Angled, fine cross-cut toothing, with cupped jaw, slimline, 17 mm long, double-action							
 Grasping and dissecting forceps	8393.288	•				•	83932887
Fine pyramid-shaped toothing, slimline, 19 mm, double-action, monopolar							
 Grasping and dissecting forceps	8393.289	•				•	83932897
Fine pyramid-shaped toothing, slimline, 23 mm, double-action							
 Biopsy forceps	8393.091	•			•		83930917
Pins, double-action			•		•		83940917
 Biopsy forceps	8393.092	•			•		83930927
Severing, single-action			•		•		83940927



∅ 5 mm, Insulated

Jaw Insert		Sheath Tube, ∅ 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
							
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
	Product No.	8393933	8394933	83930072	83930073	83930074	Bundle No.
 Bipolar forceps	8393.093	•			•		83930937
Pointed spoon forceps, single-action							
 Gastric band manipulator	8393253	•		•			83932537
	8394253		•	•			83942537
 Extracorporeal knot pusher	8393.345	•		•			83933457
Distal fork, double-action	8394345		•	•			83943457
 Scissors "Metzenbaum"	8393.041	•			•		83930417
Curved, double-action	8394.041		•		•		83940417
 Micro-scissors "Metzenbaum"	8393.044	•			•		83930447
Curved, double-action, ∅ 5 mm, monopolar							





Complete instruments in the above table are **recommendations only**. They are comprised of the components: jaw insert, sheath, and compatible handle. However, **any combination is possible within the same diameter**. The ERAGONaxial handle can also be selected as an addition or an alternative.







ø 5 mm, Insulated

Jaw Insert		Sheath Tube, ø 5 mm, Insulated		Handle			Complete Instrument • Jaw Insert • Sheath Tube • Handle
							
		WL 330 mm	WL 450 mm	With Lock, without HF	Without Lock, with HF, Mono	With Lock, with HF, Mono	
	Product No.	8393933	8394933	83930072	83930073	83930074	Bundle No.
Scissors  Slimline, straight, double-action	8393.045	•			•		83930457
Hook scissors  Serrated, single-action	8393.081	•			•		83930817
	8394.081		•		•		83940817
Scissors "Metzenbaum"  Curved, double-action	8393243	•			•		83932437
	8394243		•		•		83942437
Scissors  Fine-toothed, curved, double-action	8393246	•			•		83932467
	8394246		•		•		83942467
Scissors  Fine-toothed, fine tip, double-action	8393248	•			•		83932487
	8394248		•		•		83942487

Sheath Tube, ø 5 mm, Insulated			
			
	Product No.	WL 330 mm	WL 450 mm
ERAGON modular sheath tube	8393932	•	
 With distal cm scale	8394932		•









∅ 10 mm, Non-Insulated

		Handle	
 All jaw inserts and sheaths can also be combined with ERAGONaxial.			
		With Lock, without HF Connector	With Lock, without HF Connector, Not Rotatable
Extension Attachment	Product No.	83930082	83930085
	8988	•	




Jaw Insert	Sheath Tube, ∅ 10 mm, Non-Insulated		Handle	Complete Instrument	
				<ul style="list-style-type: none"> • Jaw Insert • Sheath Tube • Handle 	
	WL 330 mm	WL 450 mm	With Lock, without HF		
	Product No.	8395902	8396902	83930072	Bundle No.
Grasping and dissecting forceps "Maryland"	8395.225	•		•	83952257
	8396225		•	•	83962257
Double-action					
Dissecting forceps	8395.226	•		•	83952268
	8396226		•	•	83962267
Right-angled, double-action					
Grasping forceps with pike mouth	8395.236	•		•	83952368
	8396236		•	•	83962367
Double-action					

Complete instruments in the above table are **recommendations only**. They are comprised of the components: jaw insert, sheath, and compatible handle. However, **any combination is possible within the same diameter**. The ERAGONaxial handle can also be selected as an addition or an alternative.

ø 10 mm, Non-Insulated

Jaw Insert		Sheath Tube, ø 10 mm, Non-Insulated		Handle	Complete Instrument • Jaw Insert • Sheath Tube • Handle
					
		WL 330 mm	WL 450 mm	With Lock, without HF	
	Product No.	8395902	8396902	83930072	Bundle No.
Claw grasping forceps 	8395.489	•		•	83954897
Double-action					
Claw grasping forceps 	8395490	•		•	83954907
Double-action, with additional row of teeth					
Atraumatic grasping forceps "Babcock" 	8395.187	•		•	83951877
Large distal grasping surface with fine cross-cut toothings, double-action	8396187		•	•	83961877
Atraumatic grasping forceps "Babcock" 	8395.214	•		•	83952148
Double-action					
Grasping forceps 	8395.204	•		•	83952048
2 x 3 claw-shaped	8396204		•	•	83962047

∅ 10 mm, Non-Insulated

Jaw Insert		Sheath Tube, ∅ 10 mm, Non-Insulated			Handle		Complete Instrument • Jaw Insert • Sheath Tube • Handle
		WL 330 mm	WL 370 mm	WL 450 mm	With Lock, without HF	With Lock, without HF	
							
Product No.		8395902	8397902	8396902	83930072	83930082	Bundle No.
Atraumatic plate-shaped grasping forceps							
8395.234		•			•		83952347
Double-action				•	•		83962347
Spoon forceps							
8395.244		•			•		83952448
Anvil grasper							
8397.655			•			•	83976558
30° angled, double-action			•		•		83976572
Manipulator for Gastric Banding							
8397.654			•		•		83976547
90° articulating, atraumatic dissector, branch with distal slot and hook							

Complete instruments in the above table are **recommendations only**. They are comprised of the components: jaw insert, sheath, and compatible handle. However, **any combination is possible within the same diameter**. The ERAGONaxial handle can also be selected as an addition or an alternative.



Reprocessing Trays and Accessories for ERAGONmodular

Laparoscopy trays are indicated for packaging and sterilizing the following:

- Up to 13 single-piece or modular laparoscopy instruments
- Fits 3.5 mm, 5 mm, or 10 mm diameter instruments
- 1 or 2 laparoscopes and a light cable
- Up to 6 modular handles for modular instruments
- Options for miscellaneous electro-surgical tips, electrical cord (monopolar or bipolar)

Trays and Accessories

General Laparoscopy Tray , aluminum based and Radel plastic lid, with rack insert; holds forceps, scopes, and light cable 19.7" L x 8.7" W x 4.4" H.....	8585010
Long Laparoscopy Tray , aluminum based and Radel plastic lid, with rack insert; holds forceps, scopes, and light cable 24.4" L x 8.7" W x 4.4" H.....	8585012
Cleaning Brush , Ø 4 mm, TL 700 mm, disposable, 10/pkg (For 3.5 mm sheaths)	7970704
Cleaning Brush , Ø 7 mm, TL 700 mm, disposable, 10/pkg (For 5 mm sheaths)	7970707
Cleaning Brush , Ø 11 mm, TL 700 mm, disposable, 10/pkg (For 10 mm sheaths)	7970711
Water Pistol Cleaning System with 8 attachments and storage rack	6199.00
HF Cable , monopolar, 3 m, reusable, (Bovie/Valleylab/Erbe).....	8106.033
HF Unipolar Cable , Wolf to Bovie, TL 10 ft, disposable, 10/pkg.....	QD1001
Cleaning Brush , disposable, 3/pkg	40999012