

A multidisciplinary bipolar solution.



ERAGON bipolar

Bipolar Instrument System



The Bipolar Instrument System—modular, ergonomic, and safe.

As part of the ERAGON modular family, the ERAGON bipolar Instrument System is ergonomic, safe, and easy to use. These modular instruments offer a variety of jaw patterns, making them ideal for a range of multidisciplinary procedures.

The ergonomic advantage.

The ERAGON bipolar handle eases hand fatigue and enhances the surgeon's comfort with wide contact surfaces at sensitive points, simple operation of the rotation dial, and precise grip and manipulation.

Safety first.










The ERAGON bipolar Instrument System was developed with safety as a top priority. The proximal end of the jaws is positioned within the sheath to prevent tissue damage from the joint mechanism. In addition, the sheath tube provides a protected environment for the conduction of the current, eliminating the need for external insulation.

Easy to use.



The "Click-it" lock technology of the ERAGON bipolar Instrument System enables quick and easy assembly and disassembly. These instruments are also validated to be sterilized while assembled, which eases the workload and increases efficiency in the OR.





Jaw Insert		Sheath, ø 5.5 mm	Handle	Complete Instrument • Jaw Insert • Sheath • Handle
				
		WL 330 mm	Rotatable	
	Product No.	83930921	839300102	Bundle No.
Grasping and Dissecting Forceps				
	83930031	•	•	839300312
Grasping and Dissecting Forceps				
 Maryland Dissector	83930032	•	•	839300322
Grasping and Dissecting Forceps				
 Dolphin	83930033	•	•	839300332
Grasping Forceps				
 Fenestrated	83930034	•	•	839300342
Dissecting Forceps				
 Hasskamp	83930035	•	•	839300352
Metzenbaum Scissors				
	83930041	•	•	839300412

HF Bipolar Connecting Cables

Connector for Unit	Product No.	Length	Connector for Unit	Product No.	Length
2 pin Valleylab/BOWA (Int.) 	8108.133	3 m	WOLF/Martin/Berchtold/ Aesculap 	8108.131	3 m