

spirit of excellence



PRINCESS
Slender 7 mm resectoscope



The Princess is a uniquely petite resectoscope with an outer diameter of only 7 mm. Its small dimensions allow for significantly reduced dilatation, low resistance upon insertion, and precise dissection — making it ideal for use in patients with cervical stenosis or fertility issues.

## **Versatile Design for Expanded Treatment Options**

- Hybrid technology for use with bipolar or monopolar energy
- Fixed and rotatable sheath options
- Working bridge for use with instruments up to 7 Fr.
- 8 mm outer diameter available for use with a larger loop
- Compatible with the 2.7 mm, 12° or 30° PANOVIEW PLUS telescopes




<sup>\*</sup>Available with rotatable model.

<sup>\*\*</sup>BiVAP electrode can only be used with bipolar energy.



spirit of excellence





Classic resectoscope



7 mm Princess Resectoscope

## Slender 7 mm diameter for

- Reduced dilatation
- Less insertion resistance
- Decreased trauma to the cervical canal





spirit of excellence

						30	
Diameter	Model	Outer Sheath	Inner Sheath	Rotatable Inner S	Sheath Consisting of:	Working Element with Open Handle	Ø 2.7 mm, 12° Telescope
	Rigid	8653.021	8653.018				8974.412
7 mm	Rotatable	8653022	8653019	Inner sheath	86530191		
				Rotatable connector	86533643	8653.224	
	Rotatable	8655.334	8655.344	Inner sheath	8655.3441		
8 mm				Rotatable connector	8654.3642		

Electrodes	Description	Size		Product Number
		7 mm	8 mm	
Cutting electrode	Bipolar/Monopolar, 3/box, packed individually, sterile	•		4653.1313
Cutting electrode	Bipolar/Monopolar, 3/box, packed individually, sterile		•	4653.1323
Roller electrode	Bipolar/Monopolar, 3/box, packed individually, sterile	•	•	4653.023
Hook electrode	Bipolar/Monopolar, 3/box, packed individually, sterile	•	•	4653033
	Bipolar, 3/box, packed individually,			46530213
BiVAP electrode	sterile		•	46530243

## **Accessories**

Working insert for auxiliary instruments up to 7 Fr865326	84
<b>2.7 mm, 30° Telescope</b> 8974.4	02
Viewing obturator for 21 Fr8653	19
Obturator for 21 Fr	18
Obturator for 24 Fr	16
HF bipolar connection cable for Erbe VIO® 3, 3 m21196-1	19
HF monopolar connection cable for Valley lab, Bovie, Conmed, 3 m	33
Fiber light cable set Ø 2.5 mm (green color/marking green), 3 m8066253	01
Scope Protector 5/pkg8570.93	25
Universal Reprocessing Tray 22.8" L x 11.1" W x 3.8" H8585.0	23

To arrange a trial, contact your local Richard Wolf representative or call Customer Service at (800) 323-9653.