

Reprocessing simplified.



## RIWO Multifunction Trays

Universal, multi-use tray system.

# RIWO Multifunction Trays

Universal, multi-use tray system simplifies reprocessing.

RIWO Multifunction Trays are ideal for pre-cleaning, disinfection, decontamination, and rinsing of endoscopes and other instruments. The three-part tray system features a non-porous outer bin, a semi-flexible inner bin with grooves that fit snugly into the outer bin, and a transparent lid that allows you to quickly identify the contents. A small footprint makes all of the RIWO Multifunction Trays excellent choices for office and hospital use.

## Product features:

- Acting as a reservoir for cleaning and high-level disinfectants, the easy-to-clean outer bin is constructed of a non-porous material that resists discoloration through continued use.
- The semi-flexible inner bin fits securely inside the outer bin during disinfection, then is placed on top of the outer bin to drain through the bottom perforations after disinfection.
- The inner bin is easily removed for the instruments to be rinsed.
- The sturdy, but flexible, lid provides a tight seal to prevent HLD chemical fumes from escaping.
- All four RIWO Multifunction Trays have space-saving footprints that fit well in an office practice.



Instruments not included.



### Lid:

Flexible, but sturdy, construction.

Easy to remove.

Transparent to allow quick identification of instruments.

Secure seal prevents fumes from escaping.

### Inner bin:

Fits on top of outer bin to discharge solution.

Bottom perforations provide optimum drainage.

Grooves allow snug fit inside outer bin.

Easily removes to rinse instruments.



### Outer bin:

Non-porous construction resists discoloration.

Easy to clean.

Small footprint is ideal for office practices.


Instruments not included.

# RIWO Multifunction Trays

## For Disinfection

Tray	Type of Disinfection	Inner Dimensions (W x H x D)	Outer Dimensions (W x H x D)	Model Number
<p><b>RIWO Multifunction Tray for Disinfection:</b> Includes outer tray, inner tray, and lid.</p>  <p>Instruments not included.</p>	HLD Disinfection Only	432 x 100 x 150 mm	552 x 165 x 200 mm	6860.901
<p><b>RIWO Multifunction Tray for Disinfection:</b> Includes outer tray, inner tray, and lid.</p> 	HLD Disinfection Only	760 x 100 x 150 mm	881 x 165 x 200 mm	6860.911
<p><b>RIWO Multifunction Tray for Disinfection and Transport:</b> Includes tray and lid.</p> 	HLD Disinfection Only	548 x 100 x 348 mm	600 x 140 x 400 mm	509.81

## For Sterilization

Tray	Type of Disinfection	Inner Dimensions (W x H x D)	Outer Dimensions (W x H x D)	Model Number
<p><b>RIWO Multifunction Tray for Sterilization:</b> Includes outer tray, inner tray, and lid.</p>  <p>Instruments not included.</p>	Autoclavable	432 x 100 x 150 mm	552 x 165 x 200 mm	8509.01

For more information, contact your local Richard Wolf representative or call Customer Service at (800) 323-9653.