


**Custom
protection,
ultimate
security for
your valuable
flexible
endoscopes**



**Flexible
Fiberscope Tray**



Sterilize, transport and store with confidence

Soft silicone brackets secure the scope gently yet firmly. Scope easily tucks into bracket from the side rather than pushing it into a rigid bracket, risking damage.



Integrated parts box facilitates thorough sterilization and holds small components in one convenient place to avoid loss.

Stainless steel posts anchor the light cord in place.



- Avoid expensive repairs by using custom sterilization trays with brackets.
- Sturdy stainless steel construction withstands rigors of transport.
- Soft, silicone brackets safely hold the contents in place for sterilization during wrapping and storage.
- Validated for STERRAD® NX, STERRAD® 100NX, V-PRO® maX or EO.
- Field tested at one of America's premier health care facilities.

Ordering information

Item	Part Number	Dimensions
Flexible Fiberscope Tray	85850302	21.08" L x 10.39" W x 2.7" D