

Access. Diagnose. Treat.



Flexible and Semi-  
Rigid Ureteroscopes

# Comprehensive Endourological Solutions

Richard Wolf's Promise to You

A worldwide leader in endourology, Richard Wolf embraces a long-standing tradition of excellence and quality in areas such as:

- Cystoscopy
- Ureteroscopy
- Stone Therapy
- BPH
- Incontinence
- Smaller anatomy



From our rigid endoscopes to our latest chip-on-the-tip technology, we are dedicated to helping doctors improve patient outcome by providing comprehensive solutions.

# Flexible and Semi-Rigid Ureteroscopes

Access. Diagnose. Treat.



## The Needle™ 4.5 Fr. Semi-Rigid Ureteroscope

- Small 4.5 Fr. self-dilating distal tip—easy entry with less trauma to the patient
- 3 Fr. working channel allows for a variety of accessory instruments
- 315 mm and 430 mm working lengths available—suitable for adult patients and smaller anatomy




## The Ultrathin™ 6 Fr. Semi-Rigid Ureteroscope

- 6 Fr. atraumatic distal tip combines small size without sacrificing utility
- 4 Fr. working channel
- Accommodates larger biopsy and grasping accessory instruments
- Provides superior irrigation for a clearer image



## The D.O.C.™ 6.5 Fr. Dual-Channel Semi-Rigid Ureteroscope

- 6.5 Fr. atraumatic distal tip combines small size without sacrificing utility
- Dual, independent operating channels—use two instruments simultaneously
- Continuous cyclonic flow clears blood and debris for a crystal clear image



## The Viper™ Flexible Ureteroscope

- 270° up/down distal deflection permits better access to hard-to-reach areas within the upper urinary tract
- Ergonomic design and light weight offers better control and easier navigation
- Self-dilating atraumatic distal tip enables easy entry with less trauma to the patient



## The Cobra™ Dual-Channel Flexible Ureteroscope

- 270° up/down distal deflection allows for unparalleled navigation and access
- Dual 3.3 Fr. operating channels provide unmatched flow and visibility, while enabling a wide variety of therapeutic options
- Adjustable laser dial for added safety, accuracy, and ease

# Ordering Information

<b>Needle Semi-Rigid Ureteroscope</b>	
4.5/6 Fr. Needle Scope, Straight Eyepiece, 5°, 3 Fr. Insert Capacity, 315 mm WL .....	8701517
4.5/6 Fr. Needle Scope, Straight Eyepiece, 5°, 3 Fr. Insert Capacity, 430 mm WL .....	8701518
<b>Ultrathin Semi-Rigid Ureteroscope</b>	
6/7.5 Fr., 5°, 330 mm WL .....	8702.517
6/7.5 Fr., 5°, 430 mm WL .....	8702.518
<b>D.O.C. Semi-Rigid Dual-Channel Ureteroscope</b>	
6.5/8.5 Fr., Dual Channel, 5°, 330 mm WL .....	8708.517
6.5/8.5 Fr., Dual Channel, 5°, 430 mm WL .....	8708.518
<b>Accessories for D.O.C., Ultrathin, and Needle</b>	
Flexible Biopsy Forceps, 4 Fr., 600 mm .....	8734.606
Flexible Alligator Forceps, 4 Fr., 600 mm .....	8734.656
Flexible Mouse Tooth Forceps, 4 Fr., 600 mm .....	8734.686
Flexible Button Electrode, 3 Fr., 570 mm .....	823.031
Flexible Alligator Forceps, 3 Fr., 530 mm .....	828.605
Flexible Biopsy Forceps, 3 Fr., 530 mm .....	829.603
Flexible HF Button Electrode, 3 Fr., 920 mm .....	8741.80
Tray for Semi-rigid .....	8809.0015

<b>Viper Flexible Ureteroscope</b>	
6/8.8 Fr., Flexible, Ergonomic Deflection, 680 mm WL, Kit .....	7325.0761
6/8.8 Fr., Flexible, Logical Deflection, 680 mm WL, Kit .....	7325.0711
<b>Cobra Dual-Channel Flexible Ureteroscope</b>	
6/9.9 Fr., Flexible, Ergonomic Deflection, 680 mm WL, Kit .....	73260761
6/9.9 Fr., Flexible, Logical Deflection, 680 mm WL, Kit .....	73260711
<b>Accessories for Viper and Cobra</b>	
Sealing Adapter, 0-6 Fr. Capacity, Sterile, 5/box .....	20.SA65
Flexible Biopsy Forceps, 3 Fr., 920 mm .....	829.601
Flexible Grasping Forceps, 3 Fr., 920 mm .....	828.651
Flexible HF Button Electrode, 3 Fr., 920 mm .....	8741.80
Tray for Flexible Fiberscopes .....	8585.0302

Contact your local Richard Wolf representative for details.