

Targeted, Precision Instrument Cleaning



Cleaning Pistol System

Cleaning Pistol System

The Cleaning Pistol System is ideal for cleaning and drying any surgical instruments with channels or lumens. The Pistol can be used with either water or pressurized air via the pistol handle, and features these unique capabilities:






- Precision flow control
- Attachments uniquely connect to spray gun tip, designed to withstand high pressure
- Display system—ease of organization and accessibility to attachments
- Connects easily to a water faucet or plumb directly to water or air pipe
- Heat insulated grip

CLEANING PISTOL SYSTEM

6199.00	Cleaning Pistol System—packaged with 8 unique tips
---------	--

The Cleaning Pistol System complements other instrument cleaning products offered by Richard Wolf.

INSTRUMENT CLEANING PRODUCTS

Product	Catalog Number	Description and Intended Use	Quantity
	7990001	Appropriate for flexible ureteroscope, any 1-1.5 mm channel. 2 mm brush diameter, 15 mm brush length, and 1200 mm overall length	10/pkg
	7990002	Appropriate for flexible cystoscope, any 3 mm channel. 3 mm brush diameter, 20 mm brush length, and 1200 mm overall length.	10/pkg
	406041	Appropriate for an instrument with 3-4 mm channels. 3 mm brush diameter, 50 mm brush length, and 400 mm overall length	10/pkg
	468691	Appropriate for use in short channels, ports, adapters, and sealing caps. 5 mm brush diameter, 10 mm brush length, and 250 mm overall length.	10/pkg
	6162601	Pure Station™ Medical Device Pre-Cleaning System Pure Station™ with flushing (Large) 26" x 16" x 6.5" depth	1

Contact your local Richard Wolf representative to schedule an evaluation, or call Customer Service at (800) 323-9653.