

Our Flexibility is Your Security

**Perf. basket
for flexible Endoscopes**

Perf. basket for flexible Endoscopes

Thanks to its intelligent silicone holders, our perf. basket for flexible endoscopes is extremely secure and has universal application for all your flexible endoscopes from Richard Wolf.

- **For all flexible endoscopes from Richard Wolf**
- **Ideal for all low-temperature sterilization procedures** such as ethylene-oxide and plasma sterilization
- **Minimum contact surface** for optimum sterilization result
- **Protection against vibrations** provided by floating storage in silicone
- **Integrated small parts basket** for accessories
- **Safe reprocessing, storage, and transport**
- **Simple handling**



Perf. basket for flexible endoscopes

Holds: flexible endoscopes,
for sterilization (steam and low-temperature),
storage and transport,

inner dimensions (wxhxd):

432 mm x 88 mm x 272 mm,

outer dimensions (wxhxd):

448 mm x 106 mm x 288 mm33016

Fixing ring Ø 50.39mm 643301261

Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

J614en11021

