


A clean item




Reprocessing Baskets



to safely hold a telescope during transport,
storage, machine reprocessing and sterilization

Reprocessing Baskets

to safely hold a telescope during transport,
storage, machine reprocessing and sterilization



- **3 different lengths**
for virtually all endoscopic telescopes
- **Transport**
The resilient silicone support strips protect the telescope perfectly against vibration and impacts to the basket
- **Minimal support contact area**
Point support of telescope provides the conditions for the best possible cleaning



Reprocessing baskets

for machine reprocessing
and sterilization,

Width: 59 mm, Height: 54 mm,

Length:

287 mm, 255 g 38015

471 mm, 355 g 38016

610 mm, 410 g 38017

Ball fixers

for reprocessing baskets

UE = 20 pcs 64271.020

Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

J610en/230623