

Sterilization with System



RIWO-System-Tray

The universal system solution
for sterilization

RIWO-System-Tray



The universal system solution for sterilization


RIWO-System-Tray is the modular, universal system solution for transport, storage, and sterilization of your valuable instrument set.

The system offers instrument and accessory containers which can also be used individually and in combination:


A tray is created such that a lot of space is provided on two levels and is safely closed and held together with a lid.

- **Stable and long-lived**
as a result of high-quality plastic
- **For transport, storage, and sterilization**
of your complete instrument set
- **Safe and gentle fixing**
of instruments as a result of silicone peg mats
- **Storage of complete instrument sets**
in combination with the accessory container
- **Universally and flexibly**
for all specialist areas
- **Ideal for steam and low-temperature procedures**
- **Universal and safe storage**
on silicone instrument mat






	Description	Product no.
Container (WxHxD) 300 x 60 x 120 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 300 x 57 x 120 mm, outer dimensions (WxHxD): 366 x 77 x 186 mm	33019
Container (WxHxD) 300 x 60 x 200 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 300 x 57 x 200 mm, outer dimensions (WxHxD): 366 x 77 x 266 mm	33020
Container (WxHxD) 400 x 60 x 120 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 400 x 57 x 120 mm, outer dimensions (WxHxD): 466 x 77 x 186 mm	33021
Container (WxHxD) 400 x 60 x 200 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 400 x 57 x 200 mm, outer dimensions (WxHxD): 466 x 77 x 266 mm	33022
Container (WxHxD) 500 x 60 x 120 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 500 x 57 x 120 mm, outer dimensions (WxHxD): 566 x 77 x 186 mm	33023
Container (WxHxD) 500 x 60 x 200 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 500 x 57 x 200 mm, outer dimensions (WxHxD): 566 x 77 x 266 mm	33024
Container (WxHxD) 500 x 75 x 120 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 500 x 73 x 120 mm, outer dimensions (WxHxD): 566 x 92 x 186 mm	33017
Container (WxHxD) 600 x 60 x 200 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 600 x 57 x 200 mm, outer dimensions (WxHxD): 666 x 77 x 266 mm	33025
Container (WxHxD) 700 x 60 x 120 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 700 x 57 x 120 mm, outer dimensions (WxHxD): 766 x 77 x 186 mm	33026
Container (WxHxD) 700 x 60 x 200 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 700 x 57 x 200 mm, outer dimensions (WxHxD): 766 x 77 x 266 mm	33027
Container (WxHxD) 700 x 75 x 120 mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 700 x 73 x 120 mm, outer dimensions (WxHxD): 766 x 92 x 186 mm	33018



	Description	Product no.
Instrument mat for container	Instrument mat for container (corner piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 198,6 x 118,6 xmm	75219.013
Instrument mat for container	Instrument mat for container (center piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 198,6 x 158,6 xmm	75219.014
Instrument mat for container	Instrument mat for container (corner piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 118,2 x 118,2 xmm	75219.016
Instrument mat for container	Instrument mat for container (center piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 157,6 x 118,2 xmm	75219.017



	Description	Product no.
Instrument mat for lid	Instrument mat for lid, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, including 3 pushbuttons, dimensions (WxD): 172 x 62 xmm	75219.012
Instrument mat for lid	Instrument mat for lid, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, including 2 pushbuttons, dimensions (WxD): 92 x 62 xmm	75219.015
Pushbutton for lid	Pushbutton to fix lid instrument mat in place	75054.004

RIWO-System-Tray

Assignment of the containers and the corresponding instrument mats

Instrument container	Instrument mat for container	Instrument mat for lid	Pushbutton for lid
			
33019	2x 75219.016	2x 75219.015	4x 75054.004
33020	2x 75219.013	2x 75219.012	6x 75054.004
33021	2x 75219.016 1x 75219.017	3x 75219.015	6x 75054.004
33022	2x 75219.013 1x 75219.014	3x 75219.012	9x 75054.004
33023	2x 75219.016 1x 75219.017	3x 75219.015	6x 75054.004
33024	2x 75219.013 1x 75219.014	3x 75219.012	9x 75054.004
33017	2x 75219.016 1x 75219.017	3x 75219.015	6x 75054.004
33025	2x 75219.013 2x 75219.014	4x 75219.012	12x 75054.004
33026	2x 75219.016 2x 75219.017	4x 75219.015	8x 75054.004
33027	2x 75219.013 2x 75219.014	4x 75219.012	12x 75054.004
33018	2x 75219.016 2x 75219.017	4x 75219.015	8x 75054.004

