



**Hygiene ·
Instrument Reprocessing**

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - RIWO-BOX

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks





7

Rubber Seals




8

Hand blower for air insufflation




PERFORATED BASKETS

Illustration	Description	Product no.
 <p>Perforated basket for laser fiber</p> <p>also:</p>  <p>Fixing ring Ø 34.6mm</p>	<p>Holds: laser fiber, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD): 231x52x223mm, outer dimensions (WxHxD): 248x72x239mm</p> <p>for perforated basket for laser fiber, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, Ø 34,6x2,62mm, PACK=5 PCS</p>	<p>38001</p> <p>643301262</p>
 <p>Perforated basket for camera head</p> <p>also:</p> <p>Set of fixing chains</p>	<p>Holds: camera head, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD): 434x62x189mm, outer dimensions (WxHxD): 448x73x200mm</p> <p>suitable for:</p> <ul style="list-style-type: none"> 85525902: Logic HD camera head cable 3m 85525922: Logic HD camera head cable 3m 85525923: Logic HD camera head cable 5m 85525942: Logic 4K camera head 5525933: Logic HD pendul camera head cable 3m 5525833: Logic HD pendul camera head blue 5514901: Performance camera head cable 4m 5514961: Performance camera head cable 4m 5521902: Flex HD camera head cable 2,5m <p>for perforated basket for camera head, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, fixing chains: 1x 6 balls, 6x 7 balls, 1x 8 balls</p>	<p>38002</p> <p>643301174</p>
 <p>Perf. basket for video mediastinoscope</p>	<p>Holds: video mediastinoscope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with instrument holder, with small part basket, inner dimensions (WxHxD): 234x67x224mm, outer dimensions (WxHxD): 256x85x250mm</p>	<p>38007</p>

PERFORATED BASKETS




Illustration		Description	Product no.
 <p>Perf. basket for standard telescope 1</p>		<p>Holds: standard telescope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport</p> <p>inner dimensions (WxHxD): 270x52x40mm, outer dimensions (WxHxD): 297x54x59mm</p>	38015
<p>Perf. basket for standard telescope 2</p>		<p>inner dimensions (WxHxD): 454x52x40mm, outer dimensions (WxHxD): 481x54x59mm</p>	38016
<p>Perf. basket for standard telescope 3</p>		<p>inner dimensions (WxHxD): 592x52x40mm, outer dimensions (WxHxD): 619x54x59mm</p>	38017
also:	<p>Fixing chain with 4 balls</p>	<p>for perforated basket, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport</p>	64271020
	<p>Spring clip for perforated basket</p>	<p>locking device for perforated basket for telescope, exchangeable</p>	64076028
 <p>Perforated basket, upper urinary tract 1</p>		<p>Holds: ureteroscope, nephroscope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with holders for sonotrodes, with small part basket,</p> <p>inner dimensions (WxHxD): 619x65x124mm, outer dimensions (WxHxD): 635x75x132mm</p> <p>suitable for:</p> <ul style="list-style-type: none"> 8701.xxx: Ureterorenoscope 8702.xxx: Ureterorenoscope E-Line 8703.xxx: Ureterorenoscope E-Line 8704.xxx: Ureterorenoscope E-Line 8708.xxx: Continuous irrigation ureterorenoscope 8968.421: Nephroscope 12° 	38003
 <p>Perf. basket for cystoscope/telescope</p>		<p>Holds: cystoscope, telescope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with small part basket,</p> <p>inner dimensions (WxHxD): 451x65x124mm, outer dimensions (WxHxD): 467x75x132mm</p> <p>suitable for:</p> <ul style="list-style-type: none"> 8626.531: Cystoscope 12° 8632.411: Operating telescope 8643.401: PANOVIEW Plus telescope 8746.401: OP-Endoscope 12° 	38004

PERFORATED BASKETS




Illustration	Description	Product no.
 <p>Perforated basket for discoscope</p> <p>also:</p> <p>Rinsing adapter</p>	<p>Holds: discoscope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with small part basket, inner dimensions (WxHxD): 451x65x124mm, outer dimensions (WxHxD): 467x75x132mm</p> <p>suitable for:</p> <p>89210.1253: Discoscope 25° 89210.1254: Discoscope 25° 89210.3253: Discoscope 25° 89210.3254: Discoscope 25°</p>	<p>38005</p> <p>64164.039</p>
 <p>Perforated basket for hysteroscope 1</p> <p>also:</p> <p>Adapter</p>	<p>Holds: hysteroscop, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, for holding replacement sheaths, with smallpart basket, inner dimensions (WxHxD): 467x65x132mm, outer dimensions (WxHxD): 467x75x132mm</p> <p>suitable for:</p> <p>8753.411: Telescope 12°</p>	<p>38006</p> <p>64164.033</p>
 <p>Perforated basket for hysteroscope 2</p> <p>also:</p> <p>Rinsing adapter</p>	<p>Holds: Hysteroscope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with small part basket, inner dimensions (WxHxD): 451x80x124mm, outer dimensions (WxHxD): 467x90x132mm</p> <p>suitable for:</p> <p>8986.401: Hysteroscope 20° 8977.411: Mini compact hysteroscope 8793.401: PANOVIEW PLUS discoscope 8792.401: PANOVIEW PLUS discoscope 8792.411: PANOVIEW PLUS discoscope 8767.412: PANOVIEW PLUS telescope</p>	<p>38009</p> <p>64164.039</p>

X21 || Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

PERFORATED BASKETS

Illustration	Description	Product no.
 <p>Perforated basket for telescope</p> <p>also: Adapter</p>	<p>Holds: Telescope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with small part basket, inner dimensions (WxHxD): 451x65x124mm, outer dimensions (WxHxD): 467x75x132mm</p> <p>suitable for: 8920.4011: OP-endoscope bndl</p>	<p>38008</p> <p>64164.031</p>
 <p>Perforated basket for stereo telescope</p>	<p>Holds: stereo telescope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD): 451x130x189mm, outer dimensions (WxHxD): 477x139x196mm</p> <p>suitable for: 8840.402: Telescope 40° 884080: Stereo-endoscope 50°</p>	<p>38010</p>
 <p>Perforated basket for nephroscope</p> <p>also:</p> <p>Rinsing adapter</p> <p>Adapter</p>	<p>Holds: Nephroscope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with small part basket, inner dimensions (WxHxD): 451x65x164mm, outer dimensions (WxHxD): 467x75x172mm</p> <p>suitable for: 8964.401: Nephroscope 12° 8965.401: Nephroscope 20° 8965.411: Nephroscope 20° 8966.401: Nephroscope 12°</p> <p>suitable for: 8964.401/8965.401/8965.411</p> <p>suitable for: 8966.401</p>	<p>38011</p> <p>64164.038</p> <p>64164.031</p>

PERFORATED BASKETS

Illustration	Description	Product no.
 <p>Perforated basket for cysto-urethroscope</p>	<p>Holds: Cysto-urethroscope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD): 435x54x73mm, outer dimensions (WxHxD): 451x64x80mm</p> <p>suitable for: 8642.403: Cystoscope 25° 8645.403: Cystoscope 25°</p>	<p>38013</p>
 <p>Perforated basket for morcescope</p> <p>also:</p> <p>Adapter</p>	<p>Holds: Morcescope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, with small part basket, inner dimensions (WxHxD): 451x65x164mm, outer dimensions (WxHxD): 467x75x172mm</p> <p>suitable for: 8970.40: Morcescope 12° 8970402: Morcescope 12° 8970403: Morcescope 0° 8970405: Morcescope 0° 8970407: Morcescope 12</p>	<p>38012</p> <p>64164.031</p>
 <p>Perforated basket for 3D telescope</p> <p>also:</p> <p>Set of fixing chains</p>	<p>Holds: 3D telescope, for machine reprocessing, sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD): 569x69x189mm, outer dimensions (WxHxD): 583x80x200mm</p> <p>suitable for: 89346311: 3D Endoscope bndl 893446311: 3D Endoscope bndl 89346321: 3D Endoscope bndl 893446321: 3D Endoscope bndl</p> <p>for perforated basket for 3D telescope, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, fixing chains: 1x 5 balls, 1x 6 balls, 2x 7 balls, 1x 8 balls</p>	<p>38014</p> <p>643301187</p>

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - Washing - Disinfection - Neutralisation

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks

7

Sealing Caps

8

Hand Blower for air insufflation




BASKET SYSTEM

RICHARD WOLF's open mesh basket system, which has universal application, holds even the most expensive instruments safely and carefully during






- sterilization
- transport
- storage

The universal basket concept allows almost any basket configuration to be set up.













X21 || Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

BASKET SYSTEM




Illustration	Description	Product no.
 <p>Universal perforated basket 1 with instrument separator (2)</p>	<p>Holds: Instrument and accessories, for sterilization (steam and low-temperature), storage and transport, with instrument separator (2), inner dimensions (WxHxD): 709x88x259mm, outer dimensions (WxHxD): 732x117x282mm</p>	<p>33001</p>
 <p>Universal perforated basket 2 with instrument separator (4)</p>	<p>Holds: Instrument and accessories, for sterilization (steam and low-temperature), storage and transport, with instrument separator (4), inner dimensions (WxHxD): 709x186x259mm, outer dimensions (WxHxD): 743x200x282mm</p>	<p>33002</p>
 <p>Universal perforated basket 3 with instrument separator (2)</p>	<p>Holds: Instrument and accessories, for sterilization (steam and low-temperature), storage and transport, with instrument separator (2), inner dimensions (WxHxD): 513x95x233mm, outer dimensions (WxHxD): 538x117x250mm</p>	<p>33003*</p>
 <p>Universal perforated basket 4</p>	<p>Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD): 513x95x233mm, outer dimensions (WxHxD): 538x117x250mm</p>	<p>33004*</p>
 <p>Universal perforated basket 5 with instrument separator (4)</p>	<p>Holds: Instrument and accessories, for sterilization (steam and low-temperature), storage and transport, with instrument separator (4), inner dimensions (WxHxD): 513x184x233mm, outer dimensions (WxHxD): 549x200x250mm</p>	<p>33005*</p>

* suitable for standard container





BASKET SYSTEM EQUIPMENT

Illustration	Description	Product no.
 <p>Instrument mat for perforated basket</p>	Storage and transport, for holding instruments during sterilization (steam and low-temperature), dimensions (WxD): 500x226mm	33301
 <p>Instrument separator for perf. basket</p>	Holds: instruments, storage and transport, for holding instruments during sterilization (steam and low-temperature)	33302
 <p>Fastening set for perforated basket</p>	Storage and transport, for holding instruments during sterilization (steam and low-temperature)	33303
 <p>Perforated basket instrument holder</p>	Storage and transport, for aligning the axial position of instruments during sterilization (steam and low-temperature)	33304
 <p>Instrument holder for perforated basket</p>	closed 3 parts, for instruments with outer diameter, storage and transport, for holding instruments during sterilization (steam and low-temperature), outer dimensions (WxHxD): 109x40x5mm	33305
 <p>Instrument holder for perforated basket</p>	open 3 parts, storage and transport, for holding instruments during sterilization (steam and low-temperature), outer dimensions (WxHxD): 109x35x5mm	33306
 <p>Instrument holder for perforated basket</p>	closed 5 parts, storage and transport, for holding instruments during sterilization (steam and low-temperature), outer dimensions (WxHxD): 109x35x5mm	33307
 <p>Instrument holder for perforated basket</p>	open 5 parts, storage and transport, for holding instruments during sterilization (steam and low-temperature), outer dimensions (WxHxD): 109x35x5mm	33308
 <p>Instrument holder for perforated basket</p>	closed 7 parts, storage and transport, for holding instruments during sterilization (steam and low-temperature), outer dimensions (WxHxD): 109x35x5mm	33309
 <p>Instrument holder for perforated basket</p>	open 7 parts, storage and transport, for holding instruments during sterilization (steam and low-temperature), outer dimensions (WxHxD): 109x25x5mm	33310

BASKET SYSTEM EQUIPMENT


Illustration	Description	Product no.
 <p>Locking pin for mesh basket</p> <p>required for this: Knurled nut for perforated basket</p>	<p>for holding instruments during sterilization (steam and low-temperature), storage and transport</p>	<p>64019127</p> <p>64082153</p>
 <p>Set of fixing rings</p>	<p>for universal perforated basket, for holding instruments during sterilization (steam and low-temperature), storage and transport, 1x Ø 90,0x3,0mm, 1x Ø 116,0x3,0mm, PACK=2 PCS</p>	<p>74002311</p>
 <p>Perforated basket for small parts</p>	<p>Holds: Small parts, for sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD) : 109x20x109mm, outer dimensions (WxHxD) : 128x40x121mm</p>	<p>33006</p>

PERFORATED BASKETS

Illustration	Description	Product no.
 <p>Perf. basket – knee arthrosc. Basic set</p>	<p>Holds: knee arthroscopy basic set, for sterilization (steam and low-temperature), storage and transport, with instrument holders, inner dimensions (WxHxD): 513x134x233mm, outer dimensions (WxHxD): 549x150x250mm</p>	<p>33007</p>
 <p>Perf. basket – knee arthrosc. ACL/PCL 1</p>	<p>Holds: instrument set for knee arthroscopy ACL/PCL 1, for sterilization (steam and low-temperature), storage and transport, with instrument holders, inner dimensions (WxHxD): 513x134x233mm, outer dimensions (WxHxD): 549x150x250mm</p>	<p>33008</p>
 <p>Perf. basket – knee arthrosc. ACL/PCL 2</p>	<p>Holds: instrument set for knee arthroscopy ACL/PCL 2, for sterilization (steam and low-temperature), storage and transport, with instrument holders, inner dimensions (WxHxD): 513x134x233mm, outer dimensions (WxHxD): 549x150x250mm</p>	<p>33009</p>
 <p>Perf. basket For knee arthrosc. prep tab</p>	<p>Holds: knee arthroscopy preparation table, for sterilization (steam and low-temperature), storage and transport, with instrument holders, inner dimensions (WxHxD): 513x134x233mm, outer dimensions (WxHxD): 549x150x250mm</p>	<p>33010</p>




suitable for standard container

PERFORATED BASKETS

Illustration	Description	Product no.
 <p>Perf. Basket for shoulder arthroscopy</p> <p>hierzu:</p> <p>Fixing ring Ø 50.39mm</p>	<p>Holds: instrument set for shoulder arthroscopy, for sterilization (steam and low-temperature), storage and transport, with instrument holders, inner dimensions (WxHxD): 513x95x233mm, outer dimensions (WxHxD): 538x117x250mm</p> <p>for perforated basket for flexible endoscopes, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, Ø 50,39x3,53mm, PACK=10 PCS</p>	<p>33012</p> <p>643301261</p>
 <p>Perf. basket for VERTEBRIS lumbar</p>	<p>Holds: VERTEBRIS lumbar instrument set, for sterilization (steam and low-temperature), storage and transport, with instrument holders, inner dimensions (WxHxD): 513x134x233mm, outer dimensions (WxHxD): 549x150x250mm</p>	<p>33011</p>
 <p>Perf. basket for VERTEBRIS foraminotomy</p>	<p>Holds: VERTEBRIS foraminotomy instrument set, for sterilization (steam and low-temperature), storage and transport, with instrument holders,</p> <p>inner dimensions (WxHxD): 503x42x246mm, outer dimensions (WxHxD): 525x65x255mm</p>	<p>33013</p>
<p>Perforated basket for pain therapy</p>	<p>inner dimensions (WxHxD): 503x72x246mm, outer dimensions (WxHxD): 525x95x255mm</p>	<p>33014</p>

Suitable for standard container

PERFORATED BASKETS

Illustration	Description	Product no.
 <p>Perf. basket for VERTEBRIS Stenosis</p>	<p>Holds: VERTEBRIS stenosis instrument set, for sterilization (steam and low-temperature), storage and transport, with instrument holders, inner dimensions (wxhxd): 513x134x233mm, outer dimensions (wxhxd): 549x150x250mm</p>	<p>33015</p>
 <p>Perf. basket for flexible endoscopes</p> <p>hierzu:</p>  <p>Fixing ring Ø 50.39mm</p>	<p>Holds: flexible endoscopes, for sterilization (steam and low-temperature), storage and transport, inner dimensions (WxHxD): 432x88x272mm, outer dimensions (WxHxD): 448x106x288mm</p> <p>for perforated basket for flexible endoscopes, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, Ø 50,39x3,53mm, PACK=10 PCS</p>	<p>33016</p> <p>643301261</p>


||22|| Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

RIWO SYSTEM TRAYS




	Description	Product no.
Container (WxHxD) 300x60x120mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 300x57x120mm, outer dimensions (WxHxD): 366x77x186mm	33019
Container (WxHxD) 300x60x200mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 300x57x200mm, outer dimensions (WxHxD): 366x77x266mm	33020
Container (WxHxD) 400x60x120mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 400x57x120mm, outer dimensions (WxHxD): 466x77x186mm	33021
Container (WxHxD) 400x60x200mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 400x57x200mm, outer dimensions (WxHxD): 466x77x266mm	33022
Container (WxHxD) 500x60x120mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 500x57x120mm, outer dimensions (WxHxD): 566x77x186mm	33023
Container (WxHxD) 500x60x200mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 500x57x200mm, outer dimensions (WxHxD): 566x77x266mm	33024
Container (WxHxD) 500x75x120mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 500x73x120mm, outer dimensions (WxHxD): 566x92x186mm	33017
Container (WxHxD) 600x60x200mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 600x57x200mm, outer dimensions (WxHxD): 666x77x266mm	33025
Container (WxHxD) 700x60x120mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 700x57x120mm, outer dimensions (WxHxD): 766x77x186mm	33026
Container (WxHxD) 700x60x200mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 700x57x200mm, outer dimensions (WxHxD): 766x77x266mm	33027
Container (WxHxD) 700x75x120mm	Holds: Instruments, for sterilization (steam and low-temperature), storage and transport, with instrument mat, inner dimensions (WxHxD): 700x73x120mm, outer dimensions (WxHxD): 766x92x186mm	33018

RIWO SYSTEM TRAYS



	Description	Product no.
Instrument mat for container	Instrument mat for container (corner piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 198,6x118,6mm	75219.013
Instrument mat for container	Instrument mat for container (center piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 198,6x158,6mm	75219.014
Instrument mat for container	Instrument mat for container (corner piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 118,2x118,2mm	75219.016
Instrument mat for container	Instrument mat for container (center piece), for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, dimensions (WxD): 157,6x118,2mm	75219.017




	Description	Product no.
Instrument mat for lid	Instrument mat for lid, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, including 3 pushbuttons, dimensions (WxD): 172x62mm	75219.012
Instrument mat for lid	Instrument mat for lid, for fixing instruments during machine reprocessing, sterilization (steam and low-temperature), storage and transport, including 2 pushbuttons, dimensions (WxD): 92x62mm	75219.015
Pushbutton for lid	Pushbutton to fix lid instrument mat in place	75054.004

RIWO SYSTEM TRAYS

ASSIGNMENT OF THE CONTAINERS AND THE CORRESPONDING INSTRUMENT MATS

Instrument container	Instrument mat for container	Instrument mat for lid	Pushbutton for lid
			
33019	2x 75219.016	2x 75219.015	4x 75054.004
33020	2x 75219.013	2x 75219.012	6x 75054.004
33021	2x 75219.016 1x 75219.017	3x 75219.015	6x 75054.004
33022	2x 75219.013 1x 75219.014	3x 75219.012	9x 75054.004
33023	2x 75219.016 1x 75219.017	3x 75219.015	6x 75054.004
33024	2x 75219.013 1x 75219.014	3x 75219.012	9x 75054.004
33017	2x 75219.016 1x 75219.017	3x 75219.015	6x 75054.004
33025	2x 75219.013 2x 75219.014	4x 75219.012	12x 75054.004
33026	2x 75219.016 2x 75219.017	4x 75219.015	8x 75054.004
33027	2x 75219.013 2x 75219.014	4x 75219.012	12x 75054.004
33018	2x 75219.016 2x 75219.017	4x 75219.015	8x 75054.004

PROTECTIVE STERILIZATION SLEEVE

Illustration	Description	Product no.
 <p>Protective sterilization sleeve</p> <p>also:</p> <p>Plug</p>	<p>Holds: Resection, vaporization electrodes, Urethrotome, max. 5 pieces, for sterilization (steam and low-temperature), storage and transport, OD 22mm, TL 328m</p>	<p>33030</p> <p>64192.019</p>

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - Washing - Disinfection - Neutralisation

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks




7

Sealing Caps

8

Hand Blower for air insufflation

RIWO TRAYS

Illustration	Description	Product no.
 <p>Tray (WxHxD) 432x100x150mm</p>	<p>Holds: instruments, for cleaning, disinfection and neutralization, inner dimensions (WxHxD): 432x100x150mm, outer dimensions (WxHxD): 550x157x200mm</p>	<p>68601</p>
 <p>Tray (WxHxD) 760x100x150mm</p>	<p>Holds: instruments, for cleaning, disinfection and neutralization, inner dimensions (WxHxD): 760x100x150mm, outer dimensions (WxHxD): 880x157x200mm</p>	<p>68602</p>
 <p>Tray (WxHxD) 548x100x348mm</p>	<p>Holds: instruments, for cleaning, disinfection and neutralization, inner dimensions (WxHxD): 548x100x348mm, outer dimensions (WxHxD): 600x139x400mm</p>	<p>68603</p>

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - Washing - Disinfection - Neutralisation

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks


7

Sealing Caps

8





Hand Blower for air insufflation

WATER JET CLEANING PISTOL



		Product no.	
Water jet cleaning pistol incl. clean pistol mounting (A-H) and holder	for connection to water faucet with R 3/4"-thread	6199.00	
Clean gun mounting (A-H)	The adapters are suitable for cleaning the following		
	A	Sheaths and inserts, syringes, injection cannulas	15515.003
	B	Sheaths, measuring and blood pipettes	15515.004
	C	Sheaths, measuring and blood pipettes	15515.005
	D	Drainage tubes	15515.006
	E	Taps, cannulas and syringes	15515.007
	F	Douche attachments	15515.009
	G	Bottles	15515.008
H	Water jet pumps (with attachments also suitable for drying)	15515.010	

CLEANING ACCESSORIES

Illustration	Description		Product no.
	O-ring tool	for easy removal of O-rings	15106.230
	Bellows pump	for blowing air through channels	103.00
	Magnifier		9500.103
	Cleaning clamp for working elements		15242024
	Cleaning clip	for machine reprocessing of suction and irrigation handle 8475.26	8475621

ADAPTER FOR WASHER-DISINFECTOR MACHINE REPROCESSING

Illustration	Description	suitable for	Product no.
	Irrigation adapter Ø 9 mm, length 64 mm	for insulated sheath tubes 8388200 / 8388201 / 8388202 / 8388203 / 8390933 / 8391933 / 8392933 / 8393933 / 8394933 / 8392935 / 8393935 / 8394935	163502
	Rinsing adapter	for uninsulated sheath tubes	163503
	Adapter	for leakage testing in OLYMPUS endoscope washer / disinfecter	163.905
	Irrigation adapter	for flexible bronchoscope with suction valve	163912
	Adapter	for uretero reoscopes for cleaning the laser fiber channel	163914
	Irrigation adapter	for motor handpiece 8564.021	163951
	Irrigation adapter	for motor handpieces 89955.0000, 89955.0003, 89955.00001, 89955.00031	163955
	Irrigation adapter	for motor handpiece 89951.0001	163956
	Connecting part, cpl.	for flexible endoscopes	15028.083
	Connecting part, cpl.	for flexible bronchoscopes	15028103

ADAPTER FOR WASHER-DISINFECTOR MACHINE REPROCESSING









Illustration	Description	suitable for	Product no.
	Attachment	für Dichtigkeitstest in E-RDG-Spülmaschinen anderer Hersteller	15092.102
	Attachment	für Siebkorb 38008/38011/38012	64164.031
	Rinsing adapter	für Siebkorb 38011	64164.038
	Rinsing adapter	für Siebkorb 38009	64164.039
	Irrigation adapter	für DuoRotate Zangen / Scheren 5 mm / 10 mm	8399300

Illustration	for flexible urethro-cystoscopes, hysteroscopes and nasopharyngo-laryngoscopes	psuitable for	Product no.
	Luer lock connector	Ø 3 mm, männlich	886.00
	Luer lock connector Ø 4mm	Ø 4 mm, männlich Länge 29,5 mm	886.01
	Connexion de tube Ø 3mm	Ø 3 mm, weiblich	891.00
	Luer lock tube part	Ø 4 mm, weiblich	15023.041

DISPOSABLE CLEANING BRUSH FOR RIGID INSTRUMENTS




for single use (PACK=10 PCS)	Brush head Ø mm	Brush length BL mm	Total length TL mm	Color	for channel Ø mm	Product no.
Cleaning brush Ø 2mm TL 400mm	2	20	400	blue	1.0 - 1.5	7970402
Cleaning brush Ø 2mm TL 700mm			700			7970702
Cleaning brush Ø 3mm TL 400mm	3	30	400	red	1.6 - 2.5	7970403
Cleaning brush Ø 3mm TL 700mm			700			7970703
Cleaning brush Ø 4mm TL 400mm	4	30	400	yellow	2.6 - 3.5	7970404
Cleaning brush Ø 4mm TL 700mm			700			7970704
Cleaning brush Ø 5mm TL 400mm	5	48	400	green	3.6 - 4.5	7970405
Cleaning brush Ø 5mm TL 700mm			700			7970705
Cleaning brush Ø 7mm TL 400mm	7	48	400	purple	4.6 - 6.5	7970407
Cleaning brush Ø 7mm TL 700mm			700			7970707
Cleaning brush Ø 9mm TL 400mm	9	48	400	orange	6.6 - 8.5	7970409
Cleaning brush Ø 9mm TL 700mm			700			7970709
Cleaning brush Ø 11mm TL 400mm	11	48	400	white	8.6 - 10.5	7970411
Cleaning brush Ø 11mm TL 700mm			700			7970711
(PACK=5 PCS)						
Cleaning brush Ø 14mm TL 400mm	14	48	400	pink	10.6 - 13.5	7970414
Cleaning brush Ø 14mm TL 700mm			700			7970714
Cleaning brush Ø 17mm TL 400mm	17	48	400	black	13.6 - 16.5	7970417
Cleaning brush Ø 17mm TL 700mm			700			7970717
Cleaning brush Ø 20mm TL 400mm	20	48	400	light blue	16.6 - 19.5	7970420
Cleaning brush Ø 20mm TL 700mm			700			7970720

DISPOSABLE CLEANING BRUSH STORAGE RACK




- Various brushes for the entire range of Richard Wolf instruments
- Color coding
- Brushes validated for the specified diameter
- No reprocessing required



		Product no.
<p>Storage rack For cleaning brushes</p>	<p>for storing and easily removing cleaning brushes, for wall mounting, Material: stainless steel, reusable</p>	<p>7970901</p>

DISPOSABLE CLEANING BRUSH

FOR STOPCOCKS AND STOPCOCK PLUGS OF RIGID AND FLEXIBLE INSTRUMENTS




for single use (PACK= 50 ST)	Brush head BK (PACK= 50PCS)		Brush length BL		Total length TL	Color	for stopcock insert / housings with	Product no.
	BK1 Ø mm	BK2 Ø mm	BL1 mm	BL2 mm	mm		Pegs	
Double conical straight cleaning brush	conical 4-7	3	35	20	180	blue	2	7980001
	conical 5-9	4	35	20	180	red	3	7980002
	conical 6-11	5	35	20	180	yellow	4	7980003

also:



		Product no.
Dispenser for Cleaning brushes	for storing and easily removing cleaning brushes, for table mounting, Material: stainless steel, reusable	7980901


DISPOSABLE CLEANING BRUSH FOR FLEXIBLE ENDOSCOPES



for single use (PACK= 10 PCS)	Ø mm	Color	Brush length BL mm	Total length TL mm	Application		Product no.
					in working and laser channel	in the flexible endoscopes	
Cleaning brush Ø 2mm TL 1200mm	2	yellow	20	1200	Ø 1.0 - 1.5 mm	7223.xxx 7325.xxx 7326.xxx 7330.xxx 7331.xxx	7990001
Cleaning brush Ø 3mm TL 1200mm	3	green	30	1200	Ø 2.0 - 2.5 mm	7250.011 7265.xxx 7268.xxx 7269.xxx 7270.xxx 7305.xxx 7308.xxx 7309.xxx 7310.xxx	7990002
Cleaning brush Ø 0,85mm TL 1200mm	0.85	blue	10	1200	Ø 0.8 mm	7356xxx	7990003
Cleaning brush Ø 5mm TL 1200mm	5	white	20	1200	Ø 2.5 - 3.5 mm	7235xxx 7236xxx 7260xxx 7315xxx	7990005

X21 || Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

CLEANING BRUSH FOR SURFACE CLEANING



	Brush length BL mm	Product no.
Cleaning brush straight, for single use (PACK=10 PCS)	25	8691



	Brush length BL mm	Product no.
Cleaning brush angled, steam-sterilizable	42	86.90

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - Washing - Disinfection - Neutralisation

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks

7

Sealing Caps

8

Hand Blower for air insufflation

STERILIZATION ACCESSORIES / DETERGENTS AND CLEANING AGENTS



163.903



163.904

Sterilization Accessories		Product no.
Leak tester For flex. endoscopes	reusable	163.903
Pressure equalization valve	reusable	163.904

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - Washing - Disinfection - Neutralisation

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks





7


Sealing Caps

8

Hand Blower for air insufflation

STOPCOCKS

Illustration	Demountable stopcock inserts	Product no.	
	<p>Stopcock for two-way stopcock</p>	<p>15036060</p>	
	Nubs	Demountable stopcock inserts	Product no.
	<ul style="list-style-type: none"> ● ● 	<p>Stopcock bndl CAP 2.5mm capacity 1.8 mm and 2.5 mm, identification: 2 nubs (PACK= 5 PCS)</p>	<p>896.0001</p>
	<ul style="list-style-type: none"> ● ● ● 	<p>Stopcock bndl CAP 3.0mm capacity 3.0 mm, identification: 3 nubs (PACK= 5 PCS)</p>	<p>896.0002</p>
	<ul style="list-style-type: none"> ● ● ● ● 	<p>Stopcock bndl CAP 4.5mm capacity 4.2 mm, identification: 4 nubs (PACK= 5 PCS)</p>	<p>896.0003</p>



The stopcock housings and the corresponding stopcock inserts are marked with indicator nubs in order to simplify assembly.

X21 || Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - Washing - Disinfection - Neutralisation

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks

7

Sealing Caps

8

Hand Blower for air insufflation

SEALING CAPS

Illustration		Capacity Ø mm	Colour	Product no.
Sealing cap small, PACK= 10 PCS 	Sealing cap without holes	-	dark orange	89.20
	Sealing cap CAP <2,4mm	< 2,4	black	89.00
	Sealing cap CAP 2,4-3,4mm	2,4 - 3,4	red	89.01
	Sealing cap CAP 3,4-5,1mm	3,4 - 5,1	blue	89.02
	Sealing cap CAP 5,1-6,0mm	5,1 - 6,0	grey	89.03
	Sealing cap CAP 6,0-7,0mm	6,0 - 7,0	brown	89.04
	Sealing cap CAP 6,0-7,0mm	7,0 - 8,0	yellow	89.05
	Sealing cap CAP 9,5-10,1mm	9,5 - 10,1	red	89.08
	Sealing cap CAP >10,1mm	> 10,1	blue	89.09
	Sealing cap CAP 13,0-13,8mm	13,0 - 13,8	red	89.10
	Sealing cap CAP 4-6mm	4,0 - 6,0	red	89.141
	Sealing cap CAP 7,5-8,5mm	7,5 - 8,5	blue	89.161
	Sealing cap CAP 7,0mm	7,0	brown	89.171
Sealing cap with cross slot 	for trocar sleeves with thread PACK= 10 PCS	-	dark orange	89.90
RIWO rubber cap PACK= 10 PCS 	RIWO rubber cap without holes	-	dark orange	88.00
	RIWO rubber cap CAP 0,5mm	0,5	red	88.005
	RIWO rubber cap CAP 0,8mm	0,8	dark orange	88.01
	RIWO rubber cap CAP 1,4mm	1,4	blue	88.02
	RIWO rubber cap CAP 2,4mm	2,4	black	88.03
	RIWO rubber cap CAP 4,0mm	4,0	green	88.04
	Rubber cap for hysteroscope long, PACK= 1 PCS	-	dark orange	18.01

HYGIENE · INSTRUMENT REPROCESSING

1

Machine Reprocessing - Sterilization - Storage - Transport

2

Sterilization - Storage - Transport

3

Manual Reprocessing - Washing - Disinfection - Neutralisation

4

Reprocessing Accessories

5

Detergents and Cleaning Agents

6

Stopcocks







7

Sealing Caps

8

Hand Blower for air insufflation

HAND BLOWER FOR SINGLE USE FOR AIR INSUFFLATION

Illustration hand blower		Applications		
		 Pharynx and Larynx / Thoracoscopy		
		 Proktoscopy / Rektoscopy		
for single use		Ø mm	Length mm	Product no.
	Hand blower with filter for single use, sterile (PACK= 5 PCS)	49	GL 567	412700
	Rectoscopy tube Ø 21mm WL 250mm with obturator, graduated, sterile (PACK=10 PCS)	21	NL 250	4844255
	Proctoscopy tube Ø 21 WL 117mm with obturator, graduated, sterile (PACK=10 PCS)	21	NL 117	4844120

X21 || Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

