

Sterile and efficient



## **Sterile Hand Bulb**

for manual air insufflation

# Sterile Hand Bulb

for manual air insufflation

The sterile hand bulb is used as an accessory for endoscopically controlled interventions.

The bulb can be used for insufflation of sterile air to dilate the rectum during proctorectoscopy.

The device is also used in thoracoscopy to atomize talcum powder with sterile air and for the delivery of local anesthetics in the upper airways.

- Applications:
  - Proctorectoscopy
  - Pharynx and Larynx
  - Thoracoscopy
- Distal filtration attachment with membrane filter
- Disposable – simple for no complex processing
- Patient-safe application



## Hand bulb, sterile (pack of 5)

Hand bulb with filter for disposable use

Ø 49 mm, OL 567 mm .....412700

Specifications subject to change without notice.

J 613.VI.18.en.4