

Cleaning
just got easier.



**Disposable
Cleaning Brushes**

A colorful way to clean.



- Various brush sizes for use on a wide range of devices
- Color coded for ease of use
- Brushes validated for each diameter according to GLP and ISO 17025
- Disposable, to reduce risk of cross-contamination

Experience our color coded, disposable cleaning brushes, designed to offer a high level of cleaning and safety.

The sophisticated design allows for a smooth and effective cleaning of ports and long working channels—saving you time in daily clinical practice and documentation efforts.



Storage Rack

Store and easily remove cleaning brushes from wall mounted hanging system
Stainless steel, reusable.....7970901



Single-Use Surface Cleaning Brush


Brush head Ø 5 mm, brush length 42 mm and total length 159 mm, single-use, PU = 3 pcs.....40999012

Single-Use Cleaning Brushes	Brush head Ø mm	Brush head length mm	Total length mm	Color	For channel Ø mm / Fr.		Part No.
PU = 10 pcs	2	20	400	blue	1.0 - 1.5 mm	3.0 - 4.5 Fr.	7970402
			700				7970702
	3	30	400	red	1.6 - 2.5 mm	4.8 - 7.5 Fr.	7970403
			700				7970703
	4	30	400	yellow	2.6 - 3.5 mm	7.8 - 10.5 Fr.	7970404
			700				7970704
	5	48	400	green	3.6 - 4.5 mm	10.8 - 13.5 Fr.	7970405
			700				7970705
	7	48	400	purple	4.6 - 6.5 mm	13.8 - 19.5 Fr.	7970407
			700				7970707
	9	48	400	orange	6.6 - 8.5 mm	19.8 - 25.5 Fr.	7970409
			700				7970709
	11	48	400	white	8.6 - 10.5 mm	25.8 - 31.5 Fr.	7970411
			700				7970711
PU = 5 pcs	14	48	400	pink	10.6 - 13.5 mm	31.8 - 40.5 Fr.	7970414
			700				7970714
	17	48	400	black	13.6 - 16.5 mm	40.8 - 49.5 Fr.	7970417
			700				7970717
	20	48	400	light blue	16.6 - 19.5 mm	49.8 - 58.5 Fr.	7970420
			700				7970720
PU = 10 pcs	0.85	10	1200	blue	0.8 mm	2.4 Fr.	7990003
	2	15	1200	yellow	1.0 - 1.5 mm	3.0 - 4.5 Fr.	7990001
	3	20	1200	green	2.0 - 2.5 mm	6.0 - 7.5 Fr.	7990002
	5	20	1200	white	2.5 - 3.5 mm	7.5 - 10.5 Fr.	7990005

Brushes for hard to reach areas.

The color coded, double conical cleaning brushes are intended for cleaning small, hard to reach areas. The unique brush head design allows for a smooth and gentle, yet effective way to clean those small spaces.

- Easily clean stopcocks and housings
- 3 single-use brush sizes
- 2 brush heads for effective cleaning



Double conical straight cleaning brush (blue)

for stopcock inserts / housings with 2 pegs,
conical brush head: Ø 4-7 mm,
round brush head: Ø 3 mm,
PU=50 pcs, for single-use 7980001

Double conical straight cleaning brush (red)

for stopcock inserts / housings with 3 pegs,
conical brush head: Ø 5-9 mm,
round brush head: Ø 4 mm,
PU=50 pcs, for single-use 7980002

Double conical straight cleaning brush (yellow)

for stopcock inserts / housings with 4 pegs,
conical brush head: Ø 6-11 mm,
round brush head: Ø 5 mm,
PU=50 pcs, for single-use 7980003



Dispenser for cleaning brushes

for storing and easily removing cleaning brushes,
for table mounting; material: stainless steel, reusable 7980901