

Reprocessing simplified.



RIWO Multifunction Trays

Universal, multi-use tray system.

Transport. Disinfect. Store.

RIWO Multifunction Trays are ideal for transport, disinfection, and storage of medical devices and accessories. The three-part tray system features a nonporous outer bin, a semi-flexible inner bin with grooves that fit firmly into the outer bin, and a transparent lid that allows you to quickly identify the contents. A small footprint makes all of the RIWO Multifunction Trays excellent choices for office or hospital use.

- Each tray is made up of a nonporous, easy to clean material that resists discoloration through continued use during the HLD cleaning process.
- Three convenient options for storing medical devices and accessories.
- The sturdy, transparent lid allows for easy identification of contents, and provides a tight seal to prevent HLD chemical fumes from escaping.

68603


- Perfect for transportation, HLD and storage of long flexible scopes and accessories.

68602

- Uniquely designed to easily fit flexible scopes, or large quantities of rigid and semi-rigid devices and accessories.
- Perforated inner tray with grooves, providing a snug fit into the outer tray and seamless transport and cleaning




68601

- The smallest of the three trays - ideal for individual case loads of rigid and semi-rigid instruments.



*instruments not included with the purchase of a multifunction tray.

Three Multifunction Sizes

Tray		Type of Disinfection	Inner Dimensions (W x H x D)	Outer Dimensions (W x H x D)	Product No.
RIWO Multifunction Tray for Transport, Disinfection and Storage: Includes outer tray, inner tray and lid.	 <p>Instruments not included.</p>	HLD Only	432 x 100 x 150 mm	550 x 157 x 200	68601
RIWO Multifunction Tray for Transport, Disinfection and Storage: Includes outer tray, inner tray, and lid.		HLD Only	760 x 100 x 150 mm	880 x 157 x 200	68602
RIWO Multifunction Tray for Transport, Disinfection and Storage: Includes tray and lid.		HLD Only	548 x 100 x 348 mm	600 x 139 x 400	68603

For more information, contact your local Richard Wolf representative, or call Customer Service at (800) 323-9653.

Specifications subject to change without notice.

1909-09-02-1223USA