

360° rotatable, flushable
and easy to assemble



HySafe






Micro-Modular Hand Instrument System

5 Fr. Micro-Modular Hand Instrument System

The HySafe, 5 Fr. Micro-Modular Hand Instrument System features:

- A 2-piece, modular design with an integrated axial flush port to facilitate more thorough reprocessing of the entire instrument as compared to single piece instruments
- The ability to replace individual components instead of the entire instrument, providing a more cost-effective instrument solution
- Click-together technology, allowing for easy instrument assembly—in any position
- Built-in rotation knob for easy, 360° maneuverability of jaw position
- A handle with built-in overload protection to help extend the lifespan of the instrument



Jaw Pattern	Description	Outer Diameter Ø	Working Length (WL)	Product No. includes: -Insert with sheath	Handle with overload protection	Bundle No. includes: -Insert with sheath -Handle with overload protection
	Grasping Forceps (Alligator)	5 Fr.	340 mm	8642.6508		8642.6502
	Grasping Forceps (Spoon Biopsy)	5 Fr.	340 mm	8642.6008		8642.6002
	Scissors	5 Fr.	340 mm	8642.6818	8083.001	8642.6812
	Biopsy Forceps (Punch Biopsy)	5 Fr.	340 mm	8642.6318		8642.6312