

spirit of excellence

Compact, customizable in design and function





CONDOR Multi-Use Equipment Cart



The multi-functional, modular endoscopic cart solution

The Condor Multi-Use Equipment Cart offers compact storage and transportation for Richard Wolf's PIRANHA Morcellation System, including a specially designed double canister holder for enhanced operational efficiency.

- Adjustable IV pole
- Transportation handle
- Adjustable shelf (x1)
- Adjustable storage drawer
- Double canister holder
- Internal and external cable management
- Integrated medical-grade power strip with four outlets and 10 ft. power cord
- Four lockable casters
- Cart Dimensions (W x D x H): 25.1" x 23.0" x 42.6"
- Cart Weight: 120 lbs

Ordering Information

Condor Multi-Use Equipment Cart...... 2208030







spirit of excellence

The multi-functional, modular endoscopic cart solution

The Condor Multi-Use Equipment Cart enables portability and housing of Richard Wolf's endoscopic imaging systems with customizable versatility to satisfy the needs of the busy office and OR.

- Adjustable IV pole
- VESA monitor mount
- Transportation handle
- Adjustable shelf (x2)
- Adjustable storage drawer
- Internal and external cable management
- Integrated medical-grade power strip with four outlets and 10 ft. power cord
- Four lockable casters
- Cart Dimensions (W x D x H)*: 25.1" x 23.0" x 60.5"
- Cart Weight: 130 lbs

Ordering Information

Condor Multi-Use Equipment Cart.....2208030

Condor Cart Video Kit......2208029

Comprising: x1 Shelf, x1 VESA Monitor Mount,

x1 Hardware Kit



^{*}To top of monitor mount



spirit of excellence