


Custom Protection,
Ultimate Security




Scope Protectors

Scope Protectors

Custom Protection, Ultimate Security
for Your Scope Investment



The *only* scope protector validated for use during sterilization.



Validated for *all* sterilization modalities

- EO
- Prevac Steam
- All STERRAD® modalities
- All V-PRO® modalities



Versatile functionality

- Easy to apply—slide onto the scope and “click” onto the light post
- Variety of uses—sterile packaging (peel pouch or tray), unsterile storage or transport

Fits up to 5 mm scope diameter, available in two lengths

- Fits rigid scopes up to 5 mm in diameter with a straight eyepiece
- 20 cm length—ideal for arthroscopes, sinusscopes, pediatric cystoscopes
- 32 cm length—compatible with cystoscopes, hysteroscopes, laparoscopes

Ordering Information

Catalog Number	Description
8570.9205	Scope Protector, Fits up to 5 mm Straight-Eyepiece Endoscopes up to 20 cm WL, Package of 5. 
8570.9325	Scope Protector, Fits up to 5 mm Straight-Eyepiece Endoscopes up to 32 cm WL, Package of 5. 

Contact your local Richard Wolf representative to schedule an evaluation, or call Customer Service at (800) 323-9653.