

4K image capture.
Touch screen convenience.



Evolution4K Recording System

Digitally records, stores, and archives
4K images and videos

Evolution4K Recording System

4K image capture. Touch screen convenience.

Bringing digital sharpness, clarity, and convenience to the OR, **Evolution4K** Recording System is a complete digital recording and documentation system that captures up to 4K video and images from a HDMI camera source and stores captured media with patient information. It also offers the benefit of duplicating captured videos and images to USB devices for documentation, storage, and presentation.

- **Simultaneous recording** of 4K images and videos possible with one touch
- **Recording** to internal hard disk and external USB 3.0 memory possible at the same time
- **iOS wireless application** (for iPhone/iPad) for live still images, recording 4K video (via USB), surgical reports, and media review.
- **Integrated 12.7 cm touch LCD screen** provides intuitive menu navigation and quick setting changes
- **1 TB internal storage** captures up to 40 hours of 4K UHD surgical footage
- **DICOM integration:** automatic DICOM c-store (PACS storage) and c-find (MWL)*
- **Work list function** saves time by allowing cases to be added before they start

*Connectivity is dependent on existing IT infrastructure.

Technical Specifications	
Camera/Video Input Resolutions	4K & 4K UHD: 4096 x 2160 p & 3840 x 2160 p, 50/60 Hz HD: 1080 p, 1080 i, 720 p
Video File Format	AVCHD H.264 encoding .mp4 file format
Still File Format	JPEG, PNG
Video Inputs	1 x HDMI
Video Outputs	1 x HDMI
Internal Hard Drive Capacity	1 TB
Dimensions (W x H x D)	210 x 107 x 247 mm
Weight	2.72 kg (6 lbs.)

Ordering Information:

Evolution4K Recording System **EVO4K130**



Specifications subject to change without notice.

A680en150720