

## The Future of Diagnostic Thoracoscopy



### **Mini Thoracoscopy Set**

for Single Port Technique

# Mini Thoracoscopy Set

for Single Port Technique

This all-in-one Instrument Set for Single Puncture Thoracoscopy allows minimally invasive diagnostic treatment with one small incision. With the small trocar the procedure is possible under local anesthetics and without painful stress for the patient after surgery.

- Very small Thoracoscopy Set with only 5.5 mm diameter
- Optimum image quality with 50.000 pixels
- Large working channel for 3.5 mm instruments
- Very atraumatic method for the patient - reduced pain
- Suction of fluids possible with instrument inserted
- No conventional rod lenses means lens breakage minimized
- The biopsy forceps provide high quality tissue samples
- The working instruments are very easily changed and are always in the field of view
- Ergonomic design



## OP-Endoscope bndl "Leipzig Egg"

Comprising:

5.5 mm Operating Laparoscope with 3.5 mm working channel, direction of view 0°, WL 215 mm, silicate image guide (8920.401), sealing membrane, Ø 17 mm (89.103), seal (8920.311), seal cap (15176.020) .....8920.4011

## Flexible Trocar Sleeve, Ø 5.5 mm,

Plastic sheath with thread, distal straight, WL 60 mm .....8905.331  
also:

**Trocar**, spherical .....8905.3315

## Flexible Trocar Sleeve, Ø 5.5 mm,

Plastic, WL 75 mm .....8906.051  
also:

## Trocar,

Tip tapered blunt .....8906.151

## Double spoon forceps monopolar

Ø 3.5 mm bndl, modular, spoon-shaped, WL 310 mm, insulated  
Comprising:

Handle pistol shaped mono (83930074), sheath tube insulated, Ø 3.5 mm, SL 330 mm (8391933), double spoon forceps insert, Ø 3.5 mm (8391216) .....83912167

## Hook Electrode,

monopolar, Ø 3.5 mm, WL 310 mm .....8379.452

## Coagulation Button Electrode,

monopolar, Ø 3.5 mm, WL 310 mm .....8379.462

## Probe Rod,

Ø 3.5 mm, graduated, WL 310 mm .....8379.672

## Suction Tube,

Ø 3 mm, graduated, WL 450 mm .....8380.68

## Fiber light cable bndl,

Ø 2.5 mm, TL 2.3 m  
Comprising:  
Fiber light cable (80662523), adapter projector side (8095.07), adapter endoscope side (809509) .....806625231

## Perforated basket for telescope

Holds: Telescope, for machine reprocessing, sterilization, storage and transport .....38008

