


Disposable Cleaning Brush
for rigid instruments

Einmal-Reinigungsbürste
für starre Instrumente



3

	Bürstenkopf / Brush head Ø mm	Bürstenlänge BL / Brush length BL mm	Gesamtlänge GL / Total length TL mm	Farbe Color	für Kanal / for channel Ø mm	Type / Type
Reinigungsbürste zum Einmalgebrauch (VE= 10 ST) Cleaning Brush for single use (PACK= 10 PCS)	2	20	400	blau blue	1.0 - 1.5	7970402
			700			7970702
	3	30	400	rot red	1.6 - 2.5	7970403
			700			7970703
	4	30	400	gelb yellow	2.6 - 3.5	7970404
			700			7970704
	5	48	400	grün green	3.6 - 4.5	7970405
			700			7970705
	7	48	400	violett purple	4.6 - 6.5	7970407
			700			7970707
	9	48	400	orange orange	6.6 - 8.5	7970409
			700			7970709
	11	48	400	weiß white	8.6 - 10.5	7970411
			700			7970711
Reinigungsbürste zum Einmalgebrauch (VE= 5 ST) Cleaning Brush for single use (PACK= 5 PCS)	14	48	400	pink pink	10.6 - 13.5	7970414
			700			7970714
	17	48	400	schwarz black	13.6 - 16.5	7970417
			700			7970717
	20	48	400	hellblau light blue	16.6 - 19.5	7970420
			700			7970720

Storage rack see overleaf.

Aufbereitungsleiste siehe Rückseite.