


DISPOSABLE CLEANING BRUSH

FOR STOPCOCKS AND STOPCOCK PLUGS OF RIGID AND FLEXIBLE INSTRUMENTS



for single use (PACK= 50 ST)	Brush head BK (PACK= 50PCS)		Brush length BL		Total length TL	Color	for stopcock insert / housings with	Product no.
	BK1 Ø mm	BK2 Ø mm	BL1 mm	BL2 mm	mm		Pegs	
Double conical straight cleaning brush	conical 4-7	3	35	20	180	blue	2	7980001
	conical 5-9	4	35	20	180	red	3	7980002
	conical 6-11	5	35	20	180	yellow	4	7980003

also:



		Product no.
Dispenser for Cleaning brushes	for storing and easily removing cleaning brushes, for table mounting, Material: stainless steel, reusable	7980901