


DISPOSABLE CLEANING BRUSH FOR FLEXIBLE ENDOSCOPES



for single use (PACK= 10 PCS)	Ø mm	Color	Brush length BL mm	Total length TL mm	Application		Product no.
					in working and laser channel	in the flexible endoscopes	
Cleaning brush Ø 2mm TL 1200mm	2	yellow	20	1200	Ø 1.0 - 1.5 mm	7223.xxx 7325.xxx 7326.xxx 7330.xxx 7331.xxx	7990001
Cleaning brush Ø 3mm TL 1200mm	3	green	30	1200	Ø 2.0 - 2.5 mm	7250.011 7265.xxx 7268.xxx 7269.xxx 7270.xxx 7305.xxx 7308.xxx 7309.xxx 7310.xxx	7990002
Cleaning brush Ø 0,85mm TL 1200mm	0.85	blue	10	1200	Ø 0.8 mm	7356xxx	7990003
Cleaning brush Ø 5mm TL 1200mm	5	white	20	1200	Ø 2.5 - 3.5 mm	7235xxx 7236xxx 7260xxx 7315xxx	7990005