


DISPOSABLE CLEANING BRUSH FOR RIGID INSTRUMENTS



for single use (PACK=10 PCS)	Brush head Ø mm	Brush length BL mm	Total length TL mm	Color	for channel Ø mm	Product no.
Cleaning brush Ø 2mm TL 400mm	2	20	400	blue	1.0 - 1.5	7970402
Cleaning brush Ø 2mm TL 700mm			700			7970702
Cleaning brush Ø 3mm TL 400mm	3	30	400	red	1.6 - 2.5	7970403
Cleaning brush Ø 3mm TL 700mm			700			7970703
Cleaning brush Ø 4mm TL 400mm	4	30	400	yellow	2.6 - 3.5	7970404
Cleaning brush Ø 4mm TL 700mm			700			7970704
Cleaning brush Ø 5mm TL 400mm	5	48	400	green	3.6 - 4.5	7970405
Cleaning brush Ø 5mm TL 700mm			700			7970705
Cleaning brush Ø 7mm TL 400mm	7	48	400	purple	4.6 - 6.5	7970407
Cleaning brush Ø 7mm TL 700mm			700			7970707
Cleaning brush Ø 9mm TL 400mm	9	48	400	orange	6.6 - 8.5	7970409
Cleaning brush Ø 9mm TL 700mm			700			7970709
Cleaning brush Ø 11mm TL 400mm	11	48	400	white	8.6 - 10.5	7970411
Cleaning brush Ø 11mm TL 700mm			700			7970711
(PACK=5 PCS)						
Cleaning brush Ø 14mm TL 400mm	14	48	400	pink	10.6 - 13.5	7970414
Cleaning brush Ø 14mm TL 700mm			700			7970714
Cleaning brush Ø 17mm TL 400mm	17	48	400	black	13.6 - 16.5	7970417
Cleaning brush Ø 17mm TL 700mm			700			7970717
Cleaning brush Ø 20mm TL 400mm	20	48	400	light blue	16.6 - 19.5	7970420
Cleaning brush Ø 20mm TL 700mm			700			7970720