



**Hand / Fingers / Toes**

## Hand / Fingers / Toes

5.1

---

**Arthrolution Arthroscopes Overview**

---

5.2

**Telescopes and Sheaths**

---

5.3

**Punches / Forceps / Instruments for Small Joints**

---

5.4

**ENDO-CARTRIS – Instruments for Carpal Tunnel Release**

---

**arthro**  
*lution* Arthroscope

Overview

arthro <i>lution</i>			Fingers, toes	Hand	Foot, elbow		Knee, shoulder, hip
Telescope Ø			1.9 mm	2.4 mm	2.7 mm		4.0 mm
Telescopes	30°	PANOVIEW telescope	<b>8860453</b>	<b>8871413</b>	<b>8476.433</b>	<b>8672.422</b>	<b>8880.543</b>
		PANOVIEW wide angle telescope					<b>8885.443</b>
	70°	PANOVIEW telescope					<b>8880.545</b>
		PANOVIEW wide angle telescope					<b>8885.445</b>
Sheath WL			40 mm	60 mm	77.5 mm	145 mm	130 mm
Sheaths	with 1 irrigation stopcock, fixed		<b>891220124</b>	<b>891220132</b>			
	with 2 irrigation stopcocks, fixed						<b>89121.0260</b>
	with 1 irrigation stopcock, rotatable					<b>891241141</b>	
	with 2 irrigation stopcocks, rotatable			<b>891221232</b>	<b>891221241</b>	<b>891241241</b>	<b>89121.1260</b>
	High-Flow with 1 irrigation stopcock, fixed		<b>891220128</b>				
	High-Flow with 2 irrigation stopcocks, rotatable			<b>891221235</b>			<b>89121.1263</b>
	Small-caliber with 2 irrigation stopcocks, rotatable						<b>891221255</b>
Obturator	blunt				<b>891222754</b>		<b>89121.0544</b>
	semi-blunt		<b>891221952</b>	<b>891222452</b>			<b>89121.0542</b>
	conical blunt				<b>891222750</b>	<b>891242750</b>	<b>89121.0540</b>
	sharp					<b>891242761</b>	<b>89121.0640</b>
	conical pointed		<b>891221962</b>	<b>891222462</b>	<b>891222762</b>		
Sealing part for sheath						<b>891040000</b>	

Hand  
5.1

|| Subject to modifications. Please note that some products may not be available in all countries. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

## Telescope 1.9 mm / 2.4 mm for Standard Arthroscope




**AUTOCLAVE**  
**134° C / 273° F**

Direction of view		SL	Ø	Product no.
 30°	PANOVIEW <b>telescope 30° Ø 1.9mm SL 61mm</b>	61 mm	1.9 mm	<b>8860453</b>
 30°	PANOVIEW <b>telescope Ø 2.4mm SL 98mm</b>	98 mm	2.4 mm	<b>8871413</b>

Recommended accessories		Product no.
<b>Perf. basket for standard telescope 1</b>	(L x B x H) 297 x 59 x 54 mm	<b>38015</b>
<b>Fiber light cable bndl</b> Comprising: Fiber light cable (80663523), adapter, projector to light cable (8095.07), and adapter, light cable to endoscope (809509)	Ø 3.5 mm, 2.3 m long	<b>806635231</b>

**arthro** *lution* Sheath for Arthroscope

with automatic click-lock mechanism

						
<b>Sheath for arthroscope OD 2.4 mm (telescope Ø 8860453)</b>					<b>Obturator for arthroscope OD 2.4/2.8 mm semi-blunt</b>	<b>Obturator for arthroscope OD 2.4/2.8 mm conical pointed</b>
High-Flow	1 fixed stopcock	Outer Ø	WL	Product no.		
	●	2.4 mm	40 mm	<b>891220124</b>	<b>891221952</b>	<b>891221962</b>
●	●	2.8 mm		<b>891220128</b>		

Hand  
5.2




<b>Spare parts</b>		<b>Product no.</b>
<b>Irrigation ring with 1 stopcock</b>		<b>150081761</b>
<b>Stopcock bndl cap 3.0mm (PACK= 5 PCS)</b>		<b>896.0002</b>
Also: <b>Disassembly tool</b>		<b>38310.0001</b>
<b>Mounting ring</b>		<b>15008.1704</b>
<b>Spare seal set (O-rings)</b>	Comprising: O-ring, small (PACK= 5 PCS) O-ring, large (PACK= 5 PCS)	<b>15176.142</b>
Also: <b>Disassembly tool</b>		<b>15106.230</b>

VI121 || Subject to modifications. Please note that some products may not be available in all countries. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

**arthro** *lution* Sheath for Arthroscope


with automatic click-lock mechanism





Hand  
5.2


							
<b>Sheath for arthroscope OD 3.2 / 3.5 mm (telescope Ø 8871413)</b>						<b>Obturator for arthroscope OD 3.2/3.5 mm semi-blunt</b>	<b>Obturator for arthroscope OD 3.2/3.5 mm conical pointed</b>
High-Flow	1 fixed stopcock	2 rotatable stopcocks	Outer Ø	WL	Product no.		
	●		3.2 mm	60 mm	<b>891220132</b>		
		●			<b>891221232</b>	<b>891222452</b>	<b>891222462</b>
●		●	3.5 mm		<b>891221235</b>		




Spare parts		Product no.
Irrigation ring with 1 stopcock		<b>150081761</b>
Irrigation ring with 2 stopcocks		<b>150081762</b>
Stopcock bndl cap 3.0mm (PACK= 5 PCS)		<b>896.0002</b>
Also: Disassembly tool		<b>38310.0001</b>
Mounting ring		<b>15008.1704</b>
Spare seal set (O-rings)	Comprising: O-ring, small (PACK= 5 PCS) O-ring, large (PACK= 5 PCS)	<b>15176.142</b>
Also: Disassembly tool		<b>15106.230</b>

## Miniature Punch / Forceps







		Ø	WL	Product no.
<b>Hook punch</b> Ø 2 mm		2.0 mm	65 mm	<b>8486.601</b>
<b>Miniature punch</b> Ø 2.2 mm		2.2 mm	75 mm	<b>8486.602</b>
<b>Miniature forceps</b> Ø 2 mm with double spoon jaw		2.0 mm	60 mm	<b>8486.611</b>
<b>Miniature grasping forceps</b> Ø 2 mm		2.0 mm	60 mm	<b>8486.621</b>



		Sheath Ø	WL	Product no.
<b>Miniature forceps</b> Ø 2.7mm with double spoon jaw		2.7 mm	110 mm	<b>8486.614</b>
<b>Miniature grasping forceps</b> Ø 2.7mm				<b>8486.624</b>
<b>Hook punch</b> Ø 2.7mm				<b>8486.604</b>

## Instruments for Small Joints



	∅	WL	Product no.
<b>Spoon ∅ 2.4x1mm WL 100mm</b> 	2.4 x 1 mm	100 mm	<b>8871.562</b>
<b>Curette round ∅ 1.45mm WL 100mm</b> 	1.45 mm		<b>8871.571</b>
<b>Knife with handle WL 100mm</b> 	--		<b>8871.586</b>
<b>Hook probe WL 100mm H2mm</b> H = 2.0 mm 	--		<b>8871.721</b>