

Unique hygiene!



**Disposable
cleaning brushes**

Disposable cleaning brushes



Our disposable cleaning brushes offer you a high level of safety. The sophisticated concept enables brush cleaning of all Richard Wolf instruments.

Color coding helps you quickly and reliably select the validated brush for the specified instrument diameter, ensuring that you always end up with a hygienically perfect medical product.

This not only saves you time in daily clinical practice, but also a considerable amount of documentation effort.

- Various brushes for the entire range of Richard Wolf instruments
- Color coding
- Brushes validated for the specified diameter
- No reprocessing required

Disposable cleaning brushes for rigid instruments	Brush head Ø mm	Brush length mm	Total length mm	Color	For channel Ø mm	Product no.
PU = 10 pcs	2	20	400	blue	1.0 - 1.5	7970402
			700			7970702
	3	30	400	red	1.6 - 2.5	7970403
			700			7970703
	4	30	400	yellow	2.6 - 3.5	7970404
			700			7970704
	5	48	400	green	3.6 - 4.5	7970405
			700			7970705
	7	48	400	purple	4.6 - 6.5	7970407
			700			7970707
	9	48	400	orange	6.6 - 8.5	7970409
			700			7970709
	11	48	400	white	8.6 - 10.5	7970411
			700			7970711
PU = 5 pcs	14	48	400	pink	10.6 - 13.5	7970414
			700			7970714
	17	48	400	black	13.6 - 16.5	7970417
			700			7970717
	20	48	400	light blue	16.6 - 19.5	7970420
			700			7970720

For this:



Storage rack for cleaning brushes

for storing and easily removing cleaning brushes,
for wall mounting, material: stainless steel, reusable 7970901

Disposable cleaning brushes for flexible instruments	Brush head Ø mm	Brush length mm	Total length mm	Color	For channel Ø mm	Product no.
PU = 10 pcs	0.85	10	1200	blau	0.8	7990003
	2	20	1200	gelb	1.0 - 1.5	7990001
	3	30	1200	grün	2.0 - 2.5	7990002
	5	20	1200	weiß	2.5 - 3.5	7990005

The double cleaning brushes for single use are intended in particular for efficient cleaning of Richard Wolf stopcock housings and stopcock inserts. The completely removable and maintenance-free stopcock design can be cleaned highly effectively and hygienically with the help of specific brushes. The different sizes of the stopcocks were also taken into account.


- Specifically for cleaning stopcocks from Richard Wolf
- For stopcocks with 2, 3, and 4 nubs
- For single use



Dispenser for cleaning brushes

for storing and easily removing cleaning brushes,
for table mounting, material: stainless steel,
reusable

7980901



Double conical straight cleaning brush

for stopcock inserts / housings with 2 nubs,
brush head 1: conical Ø 4-7 mm,
brush head 2: Ø 3 mm,
PU=50 pcs, color blue,
for single use

7980001

Double conical straight cleaning brush

for stopcock inserts / housings with 3 nubs,
brush head 1: conical Ø 5-9 mm,
brush head 2: Ø 4 mm,
PU=50 pcs, color red,
for single use

7980002

Double conical straight cleaning brush

for stopcock inserts / housings with 4 nubs,
brush head 1: conical Ø 6-11 mm,
brush head 2: Ø 5 mm,
PU=50 pcs, color yellow,
for single use

7980003

Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

J61 Ben260826

