

spirit of excellence

## Simply smoke-free



#### Sterile Smoke-Evacuation Filter

The simple and low cost alternative to active smoke evacuation

### **Sterile Smoke-Evacuation Filter**

# The simple and low cost alternative to active smoke evacuation



spirit of excellence

As an alternative to active smoke evacuation, our sterile smoke filter offers a simple, safe and efficient option for evacuating the smoke generated while using HF current during laparoscopic interventions. The smoke is polluted with harmful particles presenting a health hazard. These pollutants are removed and captured by the filter. The filtered and treated gas exhausts to the operating room and all the operating personnel are protected from potential harm.

- Combined active charcoal and ULPA filter retains 99.999% of all particles, viruses and bacteria
- Ready to use immediately, sterile disposable
- Easy connection to the trocar with 3 mm luer-lock
- Significantly reduced odor nuisance
- Improved visual conditions in the abdomen
- Protection for the entire operating team against fine particles, viruses and bacteria present in the smoke



#### Smoke-Evacuation Filter, sterile

Richard Wolf GmbH Pforzheimer Straße 32 75438 Knittlingen, Germany Tel.: +49 7043 35-0 Fax: +49 7043 35-4300 info@richard-wolf.com

Co-CEOs: ) Jürgen Pfab Jürgen Steinbeck Trade Register: Mannheim HRB 510031 VAT no.: DE144521586 Tax no.: 48020/00171

