

Outstanding flow
properties







Optiflow





Highflow continuous irrigation
surgical sheath

Highflow continuous irrigation surgical sheath

Single continuous irrigation sheath – compact design for hysteroscopy

- Single sheath with Click-Lock quick connection for rapid assembly and additional telescope stabilization
- In combination with R. WOLF PANOVIEW hysteroscopes from the 8974 series, 30° or 12°
- Increased flow performance compared to conventional multi-sheath systems
- Intake and drainage connections angled at 45° for easy tube installation and additional working length
- Low profile instrument port for auxiliary instruments, increasing service life and enabling sensitive work
- Atraumatic instrument tip for gentle access, especially in cases involving a stenotic cervical canal
- One-piece, compact system: ergonomic and safe to use

		Product no.
	PANOVIEW telescope 30° Ø 2.7 mm SL 310 mm WL 300 mm, TL 370 mm	8974.402
	PANOVIEW telescope 12° Ø 2.7 mm SL 310 mm WL 300 mm, TL 373 mm	8974.412
	Sheath for hysteroscope 16.5 Fr WL 219 mm, round, working channel Ø 5 Fr., Optiflow sheath, distal end beveled	8974.351
	Handle for hysteroscope Ergonomic, bayonet lock	8986.251

		Product no.
	Grasping forceps bndl Pike-mouth, double-action, Ø 5 Fr., WL 340 mm	8642.6502
	Grasping forceps bndl Double-action, Ø 10.5 Fr., WL 340 mm	8642.6002
	Biopsy forceps bndl Punch, single-action, Ø 5 Fr., WL 340 mm	8642.6312
	Scissors bndl Micro scissors, single-action, Ø 5 Fr., WL 340 mm	8642.6812

Fiber light cable bndl
Length 2.3 m,
fiber bundle Ø 2.5 mm..... 806625301

Container
(W x H x D) 400 x 60 x 200 mm
Holds: Instruments, for sterilization, storage and
transport, with instrument mat..... 33022

