

A New Convenient System



Disposable Proctorectoscope Tubes

Packed sterile – for immediate use

Disposable Proctorectoscope Tubes

A New Convenient System

Richard Wolf Disposable Proctorectoscope Tubes have a long track record of providing a clear view in proctological diagnostics for many years.

This step has now become even more efficient. Sterile packing at the factory means that the proctorectoscope tubes are now ready for use immediately. This eliminates the process of disinfection prior to application.



- Improved handling**

Each tube is packed singly and sterile. Disinfection prior to application is eliminated!

- Cost-effective solution**

You benefit from the ideal combination of disposable tube and reusable handle.

- Everything in view**

An easily readable scale makes documentation and checks easier.

- Everything under control**

The reusable handle permits easy application with optimum illumination.

Tubes and accessories

	Rectoscope tube , sterile, Ø 21 mm, working length 250 mm, pack of 10 4844255
	Proctoscope tube , sterile, Ø 21 mm, working length 117 mm, pack of 10 4844120
	Handle with annular cold light illumination..... 8838.90
	Hand bulb, sterile (pack of 5) Hand bulb with filter for disposable use Ø 49 mm, OL 567 mm 412700
	Locking cap with observation window 8838.60

Subject to modifications. Please note that some products may not be available in all countries. Please always take note of the information that is provided with the product. For more detailed information, please consult with your contact partner or directly with Richard Wolf GmbH.

B 794en120225

