


spirit of excellence



Sterile Universal Sealing Valve



spirit of excellence



The new sterile universal sealing valve from Richard Wolf is a small masterpiece of engineering. It features a high level of adaptability with simple, comfortable handling. Absolute dependability and a broad range of applications make the self-closing seal system a valuable tool for streamlined and precise working.

Perfection fits everywhere

- Versatile: Ideal for diagnostics and therapy with rigid and flexible ureterorenoscopes, nephroscopes, flexible cystoscopes, etc.
- Universal: Suitable for all 3 mm luer connectors and mounting threads on ureterorenoscopes and nephroscopes. An internal and external thread on the valve makes this possible.
- Flexible: Automatically fits endoscopic auxiliary instruments up to 6 Fr.
- Comfortable: Smooth-running, self-centering applicator. Easy threading from distal and proximal end possible, without egress of irrigation fluid.
- Unique: Packed singly, sterile and useable immediately.
- Safe: Sure-grip profile for good hold even under wet and slippery conditions.









Sterile Universal Sealing Valve



For inserting auxiliary instruments of $1-6\,\,\mathrm{Fr.}$, e.g. laser fibers, LithoClast probes, tipless baskets, guide wires, forceps, etc.,

Trade Register: Mannheim HRB 510031 VAT no.: DE144521586 Tax no.: 48020/00171

